

Product Datasheet

Breast cancer suppressor candidate 1 Antibody (OTI3D6) NBP2-02898

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-02898

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-02898



NBP2-02898

Breast cancer suppressor candidate 1 Antibody (OTI3D6)


Product Information	
Unit Size	0.1 ml
Concentration	0.69 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI3D6
Preservative	0.02% Sodium Azide
Isotype	IgG2a
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	86.3 kDa

Product Description	
Description	Novus Biologicals Mouse Breast cancer suppressor candidate 1 Antibody (OTI3D6) (NBP2-02898) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	4013
Gene Symbol	VWA5A
Species	Human
Immunogen	Full length human recombinant protein of human VWA5A (NP_055437) produced in HEK293T cell.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:2000, Flow Cytometry 1:100, Immunohistochemistry 1:150, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:150

Images


Western Blot: Breast cancer suppressor candidate 1 Antibody (OTI3D6) [NBP2-02898] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY Breast cancer suppressor candidate 1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-Breast cancer suppressor candidate 1.




Immunocytochemistry/Immunofluorescence: Breast cancer suppressor candidate 1 Antibody (OTI3D6) [NBP2-02898] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY Breast cancer suppressor candidate 1.




Immunohistochemistry-Paraffin: Breast cancer suppressor candidate 1 Antibody (OTI3D6) [NBP2-02898] - Staining of paraffin-embedded Human pancreas tissue using anti-Breast cancer suppressor candidate 1 mouse monoclonal antibody.




Flow Cytometry: Breast cancer suppressor candidate 1 Antibody (OTI3D6) [NBP2-02898] - Analysis of Jurkat cells, using anti-Breast cancer suppressor candidate 1 antibody, (Red), compared to a nonspecific negative control antibody (Blue).




Immunohistochemistry-Paraffin: Breast cancer suppressor candidate 1 Antibody (OTI3D6) [NBP2-02898] - Staining of paraffin-embedded Adenocarcinoma of Human breast tissue using anti-Breast cancer suppressor candidate 1 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: Breast cancer suppressor candidate 1 Antibody (OTI3D6) [NBP2-02898] - Staining of paraffin-embedded Carcinoma of Human lung tissue using anti-Breast cancer suppressor candidate 1 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: Breast cancer suppressor candidate 1 Antibody (OTI3D6) [NBP2-02898] - Staining of paraffin-embedded Human lymph node tissue using anti-Breast cancer suppressor candidate 1 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: Breast cancer suppressor candidate 1 Antibody (OTI3D6) [NBP2-02898] - Staining of paraffin-embedded Human lymphoma tissue using anti-Breast cancer suppressor candidate 1 mouse monoclonal antibody.



Flow Cytometry: Breast cancer suppressor candidate 1 Antibody (OTI3D6) [NBP2-02898] - Analysis of Hela cells, using anti-Breast cancer suppressor candidate 1 antibody, (Red), compared to a nonspecific negative control antibody (Blue).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-02898

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
H00004013-P01-10ug	Recombinant Human Breast cancer suppressor candidate 1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-02898

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

