

Product Datasheet

GSTA4 Antibody (OTI3F5) NBP2-03188

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-03188

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-03188



NBP2-03188

GSTA4 Antibody (OTI3F5)

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI3F5
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	25.5 kDa


Product Description	
Description	Novus Biologicals Mouse GSTA4 Antibody (OTI3F5) (NBP2-03188) is a monoclonal antibody validated for use in IHC, WB and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	2941
Gene Symbol	GSTA4
Species	Human
Immunogen	Full length human recombinant protein of human GSTA4 (NP_001503) produced in HEK293T cell.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:2000, Flow Cytometry 1:100, Immunohistochemistry 1:50-1:150, Immunohistochemistry-Paraffin




Images


Western Blot: GSTA4 Antibody (OTI3F5) [NBP2-03188] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY GSTA4 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GSTA4.




Immunohistochemistry-Paraffin: GSTA4 Antibody (OTI3F5) [NBP2-03188] - Staining of paraffin-embedded Adenocarcinoma of Human breast tissue using anti-GSTA4 mouse monoclonal antibody at 1:150



Flow Cytometry: GSTA4 Antibody (OTI3F5) [NBP2-03188] - Analysis of Hela cells, using anti-GSTA4 antibody, (Red), compared to a nonspecific negative control antibody (Blue).



Immunohistochemistry-Paraffin: GSTA4 Antibody (OTI3F5) [NBP2-03188] - Staining of paraffin-embedded Human liver tissue within the normal limits using anti-GSTA4 mouse monoclonal antibody at 1:150





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-03188

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP2-52144-0.02mg	Recombinant Human GSTA4 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-03188

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

