

Product Datasheet

DCLRE1B Antibody (OTI2B7)

NBP2-03563

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-03563

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-03563



NBP2-03563**DCLRE1B Antibody (OTI2B7)**


Product Information	
Unit Size	0.1 ml
Concentration	0.62 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI2B7
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1% BSA, 50% Glycerol
Target Molecular Weight	59.8 kDa

Product Description	
Description	Novus Biologicals Mouse DCLRE1B Antibody (OTI2B7) (NBP2-03563) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	64858
Gene Symbol	DCLRE1B
Species	Human, Rat
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-280 of human DCLRE1B(NP_073747) produced in E.coli.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:200-500, Flow Cytometry 1:100, Immunohistochemistry 1:150, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin 1:150

Images

Western Blot: DCLRE1B Antibody (OTI2B7) [NBP2-03563] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DCLRE1B (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DCLRE1B.




Immunocytochemistry/Immunofluorescence: DCLRE1B Antibody (OTI2B7) [NBP2-03563] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY DCLRE1B .




Immunohistochemistry-Paraffin: DCLRE1B Antibody (OTI2B7) [NBP2-03563] - Staining of paraffin-embedded Human endometrium tissue using anti-DCLRE1B mouse monoclonal antibody.




Flow Cytometry: DCLRE1B Antibody (OTI2B7) [NBP2-03563] - HEK293T cells transfected with either overexpression plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-DCLRE1B antibody, and then analyzed by flow cytometry.




Western Blot: DCLRE1B Antibody (OTI2B7) [NBP2-03563] - 1) 293T cell lysate. 2) 293T cell transfected to overexpress 3xFlag-DCLRE1B. WB image submitted by a verified customer review.




Immunohistochemistry-Paraffin: DCLRE1B Antibody (OTI2B7) [NBP2-03563] - Staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-DCLRE1B mouse monoclonal antibody.




Immunohistochemistry-Paraffin: DCLRE1B Antibody (OTI2B7) [NBP2-03563] - Staining of paraffin-embedded Carcinoma of Human bladder tissue using anti-DCLRE1B mouse monoclonal antibody.



Flow Cytometry: DCLRE1B Antibody (OTI2B7) [NBP2-03563] - Analysis of Hela cells, using anti-DCLRE1B antibody, (Red), compared to a nonspecific negative control antibody (Blue).



Flow Cytometry: DCLRE1B Antibody (OTI2B7) [NBP2-03563] - Analysis of Jurkat cells, using anti-DCLRE1B antibody, (Red), compared to a nonspecific negative control antibody (Blue).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-03563

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-03563

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

