

# Product Datasheet

## KIR2DL1/CD158a Antibody (MM0439-8J29) [PE] NBP2-11759PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-11759PE](http://www.novusbio.com/NBP2-11759PE)

Updated 4/3/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-11759PE](http://www.novusbio.com/reviews/destination/NBP2-11759PE)



**NBP2-11759PE**

KIR2DL1/CD158a Antibody (MM0439-8J29) [PE]

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C in the dark.                                                                    |
| <b>Clonality</b>     | Monoclonal                                                                                  |
| <b>Clone</b>         | MM0439-8J29                                                                                 |
| <b>Preservative</b>  | 0.05% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG1                                                                                        |
| <b>Conjugate</b>     | PE                                                                                          |
| <b>Purity</b>        | Protein G purified                                                                          |
| <b>Buffer</b>        | PBS                                                                                         |

**Product Description**

|                                |                                                                                                                                                                                       |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>             | This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |
| <b>Host</b>                    | Mouse                                                                                                                                                                                 |
| <b>Gene ID</b>                 | 3802                                                                                                                                                                                  |
| <b>Gene Symbol</b>             | KIR2DL1                                                                                                                                                                               |
| <b>Species</b>                 | Human                                                                                                                                                                                 |
| <b>Specificity/Sensitivity</b> | This antibody only detects KIR2DL1 there is no cross reactivity with other KIRs.                                                                                                      |
| <b>Immunogen</b>               | Human recombinant KIR2DL1/CD158a                                                                                                                                                      |

**Product Application Details**

|                              |                             |
|------------------------------|-----------------------------|
| <b>Applications</b>          | Flow Cytometry, CyTOF-ready |
| <b>Recommended Dilutions</b> | Flow Cytometry, CyTOF-ready |

**Images**

Product Image: KIR2DL1/CD158a Antibody (MM0439-8J29) [PE] [NBP2-11759PE] - Vial of PE conjugated antibody. PE has two excitation maxima, 498 nm excited by the Blue laser (488 nm) and 565 nm excited by the Yellow-Green laser (561 nm). Both result in emission at 578 nm.



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP2-11759PE**

---

|              |                                                     |
|--------------|-----------------------------------------------------|
| NBP1-97005PE | Mouse IgG1 Isotype Control (MG1) [PE]               |
| NBP2-11759UV | KIR2DL1/CD158a Antibody (MM0439-8J29) [DyLight 350] |
| 202-IL-010   | IL-2 [Unconjugated]                                 |
| 1844-KR-050  | KIR2DL1/CD158a                                      |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-11759PE](http://www.novusbio.com/reviews/submit/NBP2-11759PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

