

Product Datasheet

CD23/Fc epsilon RII Antibody (MM0164-2C23) - Azide and BSA Free NBP2-12152

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-12152

Updated 2/13/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-12152



NBP2-12152

CD23/Fc epsilon RII Antibody (MM0164-2C23) - Azide and BSA Free


| Product Information | |
|-----------------------------|---|
| Unit Size | 0.1 mg |
| Concentration | LYOPH mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | MM0164-2C23 |
| Preservative | No Preservative |
| Reconstitution Instructions | Reconstitute with sterilized PBS to a final concentration of 0.5 mg/ml. |
| Isotype | IgG1 |
| Purity | Protein G purified |
| Buffer | Lyophilized from a 0.2 um filtered solution in PBS. 0.025 mg size is provided in liquid form, PBS |

| Product Description | |
|---------------------|--|
| Host | Mouse |
| Gene ID | 2208 |
| Gene Symbol | FCER2 |
| Species | Human |
| Immunogen | Human recombinant CD23/Fc epsilon RII (extracellular domain) |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot, Flow Cytometry, Block/Neutralize, CyTOF-ready |
| Recommended Dilutions | Western Blot 1:500-1:1000, Flow Cytometry 1:10-1:1000, CyTOF-ready, Block/Neutralize |
| Application Notes | This antibody is Cytof ready. |

Images

Flow Cytometry: CD23/Fc epsilon RII Antibody (MM0164-2C23) [NBP2-12152] - Analysis using the DyLight 405 conjugate of NBP2-12152. Staining of CD23 in human PBMCs using anti-CD23 antibody. Image from verified customer review.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-12152

| | |
|------------------|--|
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB720-B | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-97005-0.5mg | Mouse IgG1 Isotype Control (MG1) |
| NBP2-12152V | CD23/Fc epsilon RII Antibody (MM0164-2C23) [DyLight 405] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-12152

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

