

Product Datasheet

HSP90 Antibody (4F3.E8) [DyLight 405] NBP2-12878V

Unit Size: 0.2 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-12878V

Updated 11/9/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-12878V



NBP2-12878V

HSP90 Antibody (4F3.E8) [DyLight 405]

| Product Information | |
|----------------------------|---|
| Unit Size | 0.2 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 4F3.E8 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a |
| Conjugate | DyLight 405 |
| Purity | Protein G purified |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|--------------------------------|---|
| Description | This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |
| Host | Mouse |
| Gene ID | 3320 |
| Gene Symbol | HSP90AA1 |
| Species | Human, Mouse, Rat, Plant |
| Reactivity Notes | Immunogen displays sequence identity for non-tested species: Wheat (Triticum spp). |
| Specificity/Sensitivity | Detects approx 90kDa. This antibody detects both a and B forms of HSP90 equally well. |
| Immunogen | Recombinant Human HSP90 purified from E.coli |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details | |
|------------------------------------|---|
| Applications | ELISA, Electron Microscopy, Immunohistochemistry, Immunohistochemistry-Paraffin |
| Recommended Dilutions | ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin, Electron Microscopy |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-12878V

| | |
|--------------|---------------------------------------|
| NBP2-12878UV | HSP90 Antibody (4F3.E8) [DyLight 350] |
| NBC1-18362 | Recombinant Human HSP90 His Protein |
| 1129-ER-050 | ErbB2/Her2 [Unconjugated] |
| NBP2-76449 | Mouse HSP90 ELISA Kit (Colorimetric) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-12878V

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

