

# Product Datasheet

## TM2D1 Recombinant Protein Antigen NBP2-13436PEP

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-13436PEP](http://www.novusbio.com/NBP2-13436PEP)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-13436PEP](http://www.novusbio.com/reviews/destination/NBP2-13436PEP)



**NBP2-13436PEP**

## TM2D1 Recombinant Protein Antigen

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Preservative	No Preservative
Purity	>80% by SDS-PAGE and Coomassie blue staining
Buffer	PBS and 1M Urea, pH 7.4.
Target Molecular Weight	25 kDa

Product Description	
Description	<p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human TM2D1.</p> <p><b>Source:</b> <i>E. coli</i></p> <p><b>Amino Acid Sequence:</b> CEDLKVGQYICKDPKINDATQEPVNCTNYTAHVSCFPAPNITCKDSSGNETHFT GNEVGFFKPISCRNVNG</p> <p><b>Fusion Tag:</b> N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G)</p> <p><b>This product is intended to be used as a blocking antigen for antibody competition assays. Any other use of this antigen is done at the risk of the user. The use of this product for commercial production is strictly prohibited. Please contact technical support if you have any questions.</b></p>
Gene ID	83941
Gene Symbol	TM2D1
Species	Human

Product Application Details	
Applications	Antibody Competition
Recommended Dilutions	Antibody Competition 10 - 100 molar excess

**Application Notes**

This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP2-13436.

It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml.

For current lot information, including availability, please contact our technical support team click [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)

For further blocking peptide related information and a protocol, click [here](#).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-13436PEP**

---

H00083941-P01-10ug	Recombinant Human TM2D1 GST (N-Term) Protein
355-BM-010	BMP-2 [Unconjugated]
NBP2-13436	TM2D1 Antibody
NB300-269	Fibrillarin Antibody (38F3)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-13436PEP](http://www.novusbio.com/reviews/submit/NBP2-13436PEP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

