

Product Datasheet

preproANP Antibody - Azide Free

NBP2-14873

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 6

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-14873

Updated 9/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-14873



NBP2-14873


preproANP Antibody - Azide Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 1% BSA, 20% Glycerol
Target Molecular Weight	17 kDa
Product Description	
Description	Novus Biologicals Rabbit preproANP Antibody - BSA Free (NBP2-14873) is a polyclonal antibody validated for use in IHC and WB. Anti-preproANP Antibody: Cited in 6 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	4878
Gene Symbol	NPPA
Species	Mouse, Rat
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: Canine (83%), Porcine (87%), Rabbit (81%), Bovine (83%).
Immunogen	Full length human preproANP Recombinant protein.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:500-1:20000




Images


Western Blot: preproANP Antibody [NBP2-14873] - Various tissue extracts (50 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with ANP antibody [N1C3] diluted at 1:2000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Western Blot: preproANP Antibody [NBP2-14873] - Knockdown of ATE1 attenuates cardiac hypertrophy. Western blot analysis of preproANP protein levels in heart samples from sham and renal ligated rats (B2) bar graph was generated by quantifying the blots and normalizing the intensities of bands to the untreated lane (C) BNP and (D) beta-MHC in Ligated ATE1 siRNA when compared with Ligated NS siRNA samples. Statistical analysis was carried out by one-way ANOVA. Data are represented as mean \pm SE. Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/s41598-019-57379-7>) licensed under a CC-BY license.



Western Blot: preproANP Antibody [NBP2-14873] - Knockdown of ATE1 attenuates cardiac hypertrophy. Western blot analysis of preproANP protein levels in heart samples from sham and renal ligated rats (B2) bar graph was generated by quantifying the blots and normalizing the intensities of bands to the untreated lane (C) BNP and (D) beta-MHC in Ligated ATE1 siRNA when compared with Ligated NS siRNA samples. Statistical analysis was carried out by one-way ANOVA. Data are represented as mean \pm SE. Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/s41598-019-57379-7>) licensed under a CC-BY license.



Publications

Li S, Cao W, Wang B et al. TRIF/miR-34a mediates aldosterone-induced cardiac inflammation and remodeling Clin. Sci. 2020-06-26 [PMID: 32542395]

Singh K, Gupta A, Sarkar A et al. Arginyltransferase knockdown attenuates cardiac hypertrophy and fibrosis through TAK1-JNK1/2 pathway Sci Rep 2020-01-17 [PMID: 31953451] (WB, Rat)

Das NA, Carpenter AJ, Yashida T et al. TRAF3IP2 mediates TWEAK/TWEAKR-induced pro-fibrotic responses in cultured cardiac fibroblasts and the heart J Mol Cell Cardiol. 2018-07-05 [PMID: 29981796] (WB, Mouse)

Zhang QJ, Tran TAT, Wang M et al. Histone lysine dimethyl-demethylase KDM3A controls pathological cardiac hypertrophy and fibrosis Nat Commun 2018-12-07 [PMID: 30531796] (WB, Mouse)

Yariswamy M, Yoshida T, Valente AJ et al. Cardiac-Restricted Overexpression of TRAF3 Interacting Protein 2 (TRAF3IP2) Results in Spontaneous Development of Myocardial Hypertrophy, Fibrosis and Dysfunction. J Biol Chem. [PMID: 27466370] (WB, Mouse)

Somanna NK, Yariswamy M, Garagliano JM et al. Aldosterone-induced cardiomyocyte growth, and fibroblast migration and proliferation are mediated by TRAF3IP2 Cell. Signal. 2015-07-04 [PMID: 26148936] (WB, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-14873

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00004878-P01-10ug	Recombinant Human preproANP GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-14873

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

