

Product Datasheet

COP9 signalosome complex subunit 2 Antibody - BSA Free NBP2-15958

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-15958

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-15958



NBP2-15958

COP9 signalosome complex subunit 2 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol
Target Molecular Weight	52 kDa

Product Description	
Description	Novus Biologicals Rabbit COP9 signalosome complex subunit 2 Antibody - BSA Free (NBP2-15958) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	9318
Gene Symbol	COPS2
Species	Human, Mouse
Immunogen	Recombinant protein encompassing a sequence within the center region of human COP9 signalosome complex subunit 2. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:5000-1:20000, Immunohistochemistry 10 - 1:500, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Western Blot: COP9 signalosome complex subunit 2 Antibody [NBP2-15958] - Sample (30 ug of whole cell lysate) A: NIH-3T3 B: JC C: BCL-1 10% SDS PAGE gel, diluted at 1:5000.




Immunocytochemistry/Immunofluorescence: COP9 signalosome complex subunit 2 Antibody [NBP2-15958] - COP9 signalosome complex subunit 2 Antibody detects COP9 signalosome complex subunit 2 protein at nucleus by immunofluorescent analysis. Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: COP9 signalosome complex subunit 2 protein stained by COP9 signalosome complex subunit 2 antibody diluted at 1:500. Blue: Hoechst 33342 staining.




Immunohistochemistry-Paraffin: COP9 signalosome complex subunit 2 Antibody [NBP2-15958] - COP9 signalosome complex subunit 2 Antibody detects COP9 signalosome complex subunit 2 protein at cytosol and nucleus on mouse colon by immunohistochemical analysis. Sample: Paraffin-embedded mouse colon. COP9 signalosome complex subunit 2 antibody dilution: 1:500.




Western Blot: COP9 signalosome complex subunit 2 Antibody [NBP2-15958] - Sample (30 ug of whole cell lysate) A: 293T B: A431 C: HeLa D: HepG2 10% SDS PAGE gel, diluted at 1:10000.



Immunohistochemistry-Paraffin: COP9 signalosome complex subunit 2 Antibody [NBP2-15958] - Immunohistochemical analysis of paraffin-embedded U87 xenograft, using antibody at 1:500 dilution.



Immunohistochemistry-Paraffin: COP9 signalosome complex subunit 2
Antibody [NBP2-15958] - COP9 signalosome complex subunit 2
Antibody detects COP9 signalosome complex subunit 2 protein at
cytosol and nucleus on mouse colon by immunohistochemical analysis.
Sample: Paraffin-embedded mouse colon. COP9 signalosome complex
subunit 2 antibody dilution: 1:500.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-15958

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00009318-P01-10ug	Recombinant Human COP9 signalosome complex subunit 2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-15958

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

