

Product Datasheet

PDIA6 Antibody NBP2-19750

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-19750

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-19750




NBP2-19750**PDIA6 Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol, 1% BSA
Target Molecular Weight	48 kDa
Product Description	
Description	Novus Biologicals Rabbit PDIA6 Antibody (NBP2-19750) is a polyclonal antibody validated for use in IHC, WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	10130
Gene Symbol	PDIA6
Species	Human, Mouse, Rat
Reactivity Notes	Chimpanzee (99%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human PDIA6. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunoprecipitation 1:100-1:500, Immunohistochemistry-Paraffin 1:100-1:1000




Images


Western Blot: PDIA6 Antibody [NBP2-19750] - Sample (30 ug of whole cell lysate) A: Raji 10% SDS PAGE diluted at 1:1000




Immunohistochemistry-Paraffin: PDIA6 Antibody [NBP2-19750] - Paraffin-embedded rat liver. PDIA6 antibody [N1N3] diluted at 1:500.



Immunohistochemistry-Paraffin: PDIA6 Antibody [NBP2-19750] - Paraffin-embedded mouse liver. PDIA6 antibody [N1N3] diluted at 1:500.



Immunoprecipitation: PDIA6 Antibody [NBP2-19750] - Immunoprecipitation of PDIA6 protein from 293T whole cell extracts using 5 ug of PDIA6 antibody [N1N3]. Western blot analysis was performed using PDIA6 antibody [N1N3]. EasyBlot anti-Rabbit IgG was used as a secondary reagent.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-19750

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-52142-0.02mg	Recombinant Human PDIA6 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-19750

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

