

Product Datasheet

UBE2O Antibody NBP2-20787

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-20787

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-20787



NBP2-20787


UBE2O Antibody

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol, 1% BSA
Target Molecular Weight	51 kDa
Product Description	
Description	Novus Biologicals Rabbit UBE2O Antibody (NBP2-20787) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	63893
Gene Symbol	UBE2O
Species	Human, Mouse
Immunogen	Recombinant protein encompassing a sequence within the center region of human UBE2O. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 10 - 1:500, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000



Images


Immunocytochemistry/Immunofluorescence: UBE2O Antibody [NBP2-20787] - UBE2O antibody [C1C2], Internal detects UBE2O protein at cytosol of methanol-fixed A431, using antibody at 1:500 dilution.




Immunohistochemistry-Paraffin: UBE2O Antibody [NBP2-20787] - UBE2O antibody [C1C2], Internal detects UBE2O protein at cytosol on mouse heart by immunohistochemical analysis. Sample: Paraffin-embedded mouse heart. UBE2O antibody [C1C2], Internal dilution: 1:500.



Immunohistochemistry-Paraffin: UBE2O Antibody [NBP2-20787] - Immunohistochemical analysis of paraffin-embedded breast ca, using antibody at 1:500 dilution.



Western Blot: UBE2O Antibody [NBP2-20787] - Various whole cell extracts (30 ug) were separated by 5% SDS-PAGE, and the membrane was blotted with UBE2O antibody [C1C2], Internal (NBP2-20787) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-20787

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00063893-P01-10ug	Recombinant Human UBE2O GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-20787

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

