

Product Datasheet

UQCRC1 Antibody - BSA Free

NBP2-20815

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-20815

Updated 2/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-20815



NBP2-20815


UQCRC1 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol
Target Molecular Weight	53 kDa
Product Description	
Host	Rabbit
Gene ID	7384
Gene Symbol	UQCRC1
Species	Human
Reactivity Notes	Mouse (88%), Rat (87%), Rhesus Monkey (97%), Chimpanzee (100%), Bovine (88%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human UQCRC1. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:1000-1:20000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000




Images


Western Blot: UQCRC1 Antibody [NBP2-20815] - A Huh7 whole cell lysate/extract B Hep3B whole cell lysate/extract C HepG2 whole cell lysate/extract 10 % SDS-PAGE UQCRC1 antibody dilution: 1:5000




Immunocytochemistry/Immunofluorescence: UQCRC1 Antibody [NBP2-20815] - HepG2 cells were fixed in 2% paraformaldehyde/culture medium at 37 (c) for 30 min. Green: NYREN18 protein stained by NYREN18 antibody [N1N3] diluted at 1:500. Blue: Hoechst 33342 staining.



Western Blot: UQCRC1 Antibody [NBP2-20815] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with UQCRC1 antibody.



Western Blot: UQCRC1 Antibody [NBP2-20815] - A. 30 ug 293T whole cell lysate/extract B. 30 ug A431 whole cell lysate/extract C. 30 ug HeLa whole cell lysate/extract D. 30 ug HepG2 whole cell lysate/extract E. 30 ug A375 whole cell lysate/extract 10 % SDS-PAGE UQCRC1 antibody dilution: 1:10000





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-20815

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-51608-0.1mg	Recombinant Human UQCRC1 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-20815

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

