

# Product Datasheet

## SNRPD3 Antibody - BSA Free

### NBP2-22305

Unit Size: 100 ul

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-22305](http://www.novusbio.com/NBP2-22305)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-22305](http://www.novusbio.com/reviews/destination/NBP2-22305)



**NBP2-22305**

SNRPD3 Antibody - BSA Free

| Product Information |                                       |
|---------------------|---------------------------------------|
| Unit Size           | 100 ul                                |
| Concentration       | 1.0 mg/ml                             |
| Storage             | Store at 4C. Do not freeze.           |
| Clonality           | Polyclonal                            |
| Preservative        | 0.09% Sodium Azide                    |
| Isotype             | IgG                                   |
| Purity              | Immunogen affinity purified           |
| Buffer              | Tris-Citrate/Phosphate (pH 7.0 - 8.0) |


| Product Description |                                                                                                                                                                                                               |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Novus Biologicals Rabbit SNRPD3 Antibody - BSA Free (NBP2-22305) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.                |
| Host                | Rabbit                                                                                                                                                                                                        |
| Gene ID             | 6634                                                                                                                                                                                                          |
| Gene Symbol         | SNRPD3                                                                                                                                                                                                        |
| Species             | Human                                                                                                                                                                                                         |
| Immunogen           | The immunogen this antibody was made to, maps to a region between residue 1 to 50 of human Small Nuclear Ribonucleoprotein D3 Polypeptide 18kDa using the numbering given in entry NP_004166.1 (GeneID 6634). |

| Product Application Details |                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Immunoprecipitation (Negative)                                               |
| Recommended Dilutions       | Western Blot 1:1000 - 1:5000, Immunohistochemistry 1:500 - 1:2000, Immunohistochemistry-Paraffin 1:500 - 1:2000, Immunoprecipitation (Negative) |
| Application Notes           | Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.                                                            |




## Images

Western Blot: SNRPD3 Antibody [NBP2-22305] - Whole cell lysate from HeLa (15 and 50 ug), 293T (50 ug), and Jurkat (50 ug) cells. NBP2-22305 used for WB at 0.4 ug/ml. Detection: Chemiluminescence with an exposure time of 3 minutes.



Immunohistochemistry-Paraffin: SNRPD3 Antibody [NBP2-22305] - Sample: FFPE section of human breast carcinoma. Antibody: Affinity purified rabbit anti- Sm-D3 used at a dilution of 1:1,000 (1ug/ml). Detection: DAB. Counterstain: hematoxylin (blue).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-22305**

---

|                 |                                                     |
|-----------------|-----------------------------------------------------|
| HAF008          | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160          | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891      | Rabbit IgG Isotype Control                          |
| NBP1-78811-50ug | Recombinant Human SNRPD3 His Protein                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-22305](http://www.novusbio.com/reviews/submit/NBP2-22305)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

