

# Product Datasheet

## GIT1 Antibody (S39B-8) [PerCP] NBP2-22423PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-22423PCP](http://www.novusbio.com/NBP2-22423PCP)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-22423PCP](http://www.novusbio.com/reviews/destination/NBP2-22423PCP)



**NBP2-22423PCP**

GIT1 Antibody (S39B-8) [PerCP]

| Product Information         |                                                                                                                                                                             |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                 |
| Storage                     | Store at 4C in the dark.                                                                                                                                                    |
| Clonality                   | Monoclonal                                                                                                                                                                  |
| Clone                       | S39B-8                                                                                                                                                                      |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                          |
| Isotype                     | IgG1                                                                                                                                                                        |
| Conjugate                   | PerCP                                                                                                                                                                       |
| Purity                      | Protein G purified                                                                                                                                                          |
| Buffer                      | PBS                                                                                                                                                                         |
| Product Description         |                                                                                                                                                                             |
| Host                        | Mouse                                                                                                                                                                       |
| Gene ID                     | 28964                                                                                                                                                                       |
| Gene Symbol                 | GIT1                                                                                                                                                                        |
| Species                     | Human, Mouse, Rat                                                                                                                                                           |
| Reactivity Notes            | Mouse reactivity reported in scientific literature (PMID: 30546041).                                                                                                        |
| Specificity/Sensitivity     | Detects approx 90kDa. Does not cross-react with GIT2.                                                                                                                       |
| Immunogen                   | Fusion protein amino acids 375-770 (C-terminus) of rat GIT1                                                                                                                 |
| Product Application Details |                                                                                                                                                                             |
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation, Proximity Ligation Assay, Knockout Validated |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Immunohistochemistry-Frozen, Proximity Ligation Assay, Knockout Validated |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-22423PCP**

---

|               |                                          |
|---------------|------------------------------------------|
| NBP1-97005PCP | Mouse IgG1 Isotype Control (MG1) [PerCP] |
| NBP1-86144PEP | GIT1 Recombinant Protein Antigen         |
| DY8485-05     | GIT1 [Biotin]                            |
| NB300-605     | iNOS Antibody - BSA Free                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-22423PCP](http://www.novusbio.com/reviews/submit/NBP2-22423PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

