

Product Datasheet

SHIP2/INPPL1 [p Tyr497] Antibody NBP2-24461

Unit Size: 0.02 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-24461

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-24461



NBP2-24461**SHIP2/INPPL1 [p Tyr497] Antibody**

Product Information	
Unit Size	0.02 mg
Concentration	0.2 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS and 0.05% BSA


Product Description	
Description	Novus Biologicals Rabbit SHIP2/INPPL1 [p Tyr497] Antibody (NBP2-24461) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3636
Gene Symbol	INPPL1
Species	Human, Mouse, Rat, Bovine, Canine, Equine, Primate
Immunogen	A portion of human SHIP2 protein containing a phosphorylated tyrosine residue at position 497 was used as the immunogen.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 0.01 - 1.0 ug/ml, Immunohistochemistry 5 ug/ml, Immunohistochemistry-Paraffin




Images

Western Blot: SHIP2/INPPL1 [p Tyr497] Antibody [NBP2-24461] - Analysis of phosphorylated SHIP2 (Y497) in recombinant fusion protein containing 1) a phosphorylated tyrosine at position 497 and 2) an unphosphorylated tyrosine at position 497, using SHIP2 phospho antibody at 0.01 ug/ml.



Immunohistochemistry-Paraffin: SHIP2/INPPL1 [p Tyr497] Antibody [NBP2-24461] - Analysis of phosphorylated SHIP2 in formalin-fixed, paraffin-embedded human brain tissue using an isotype control (top) and this antibody (bottom) at 5 ug/ml.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-24461

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00003636-Q01-10ug	Recombinant Human SHIP2/INPPL1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-24461

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

