

Product Datasheet

STING/TMEM173 Antibody

NBP2-24548

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 5

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-24548

Updated 2/8/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-24548



NBP2-24548**STING/TMEM173 Antibody**

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Purity	Immunogen affinity purified
Buffer	0.2 ml PBS and 0.05% BSA

Product Description	
Host	Rabbit
Gene ID	340061
Gene Symbol	TMEM173
Species	Human
Reactivity Notes	The amino acid sequence used as immunogen is 100% homologous in human (isoform CRA_a), chimpanzee and rhesus monkey, 87% homologous in dog and cow, 81% homologous in horse, mouse and boar and 75% homologous in rat.
Immunogen	A portion of amino acids 50-100 of human STING was used as the immunogen.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 5-10 ug/ml, Immunocytochemistry/ Immunofluorescence

Publications

Diner BA, Li T, Greco TM et al. The functional interactome of PYHIN immune regulators reveals IFIX is a sensor of viral DNA. Mol Syst Biol. 2015-02-09 [PMID: 25665578]

Berg RK, Rahbek SH, Kofod-Olsen E et al. T Cells Detect Intracellular DNA but Fail to Induce Type I IFN Responses: Implications for Restriction of HIV Replication. PLoS One 2014-01-01 [PMID: 24404168] (WB, ICC/IF, Human)

Holm CK, Jensen SB, Jakobsen MR et al. Virus-cell fusion as a trigger of innate immunity dependent on the adaptor STING. Nat Immunol. 2012-06-17 [PMID: 22706339]

Tamassia N, Bazzoni F, Le Moigne V et al. IFN- β expression is directly activated in human neutrophils transfected with plasmid DNA and is further increased via TLR-4-mediated signaling. J Immunol. 2012-08-01 [PMID: 22730532] (WB)

Details:

TBk1 (IMG-139A) & STING (IMG- 6422A) 1. WB (TBk1): Whole human neutrophil and monocytes, Fig S1. 2. WB (STING): Whole human neutrophil and monocytes, Fig 6A.

Heiber JF, Barber GN. Evaluation of innate immune signaling pathways in transformed cells. Methods Mol Biol. 2012-01-01 [PMID: 21948479] (WB)

Details:

WB (Table 1), page 226. Data not shown.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-24548

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

