

# Product Datasheet

## STING/TMEM173 Antibody [PE/Atto594] NBP2-24683PEATT594

Unit Size: 0.1 ml

Store at 4C in the dark. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-24683PEATT594](http://www.novusbio.com/NBP2-24683PEATT594)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-24683PEATT594](http://www.novusbio.com/reviews/destination/NBP2-24683PEATT594)



**NBP2-24683PEATT594**

STING/TMEM173 Antibody [PE/Atto594]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark. Do not freeze.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	PE/Atto594
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	42 kDa

Product Description	
Host	Rabbit
Gene ID	340061
Gene Symbol	STING1
Species	Human, Mouse, Canine, Primate, Rhesus Macaque
Reactivity Notes	Opossum, Zebrafish (83%), Xenopus (72%), Rat (88%).
Immunogen	Partial synthetic peptide made to an internal portion of human STING/TMEM173 (between amino acids 310-360) [UniProt Q86WV6]

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen
Application Notes	Optimal dilution of this antibody should be experimentally determined.

**Images**

Product Image: STING/TMEM173 Antibody [PE/Atto594] [NBP2-24683PEATT594] - Vial of PE/Atto 594 conjugated antibody. PE/Atto 594 has two excitation maxima, 498 nm excited by the Blue laser (488 nm) and 566 nm excited by the Yellow-Green laser (561 nm). Both result in emission at 632 nm.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-24683PEATT594**

---

NBP2-24683PEP	STING/TMEM173 Antibody Blocking Peptide
210-TA-005	TNF-alpha [Unconjugated]
NBL1-17038	STING/TMEM173 Overexpression Lysate
NBP2-24729	TLR9 Antibody (26C593.2) - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-24683PEATT594](http://www.novusbio.com/reviews/submit/NBP2-24683PEATT594)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

