

Product Datasheet

LOX Antibody - BSA Free NBP2-24877

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 23

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-24877

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-24877



NBP2-24877

LOX Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Sodium Azide
Isotype	IgG
Purity	Unpurified
Buffer	Whole sera with 50% Glycerol


Product Description	
Description	Novus Biologicals Knockout (KO) Validated Rabbit LOX Antibody - BSA Free (NBP2-24877) is a polyclonal antibody validated for use in IHC, WB, ICC/IF and IP. Anti-LOX Antibody: Cited in 22 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	4015
Gene Symbol	LOX
Species	Human, Mouse, Rat, Porcine, Sheep
Reactivity Notes	Porcine reactivity reported in scientific literature (PMID: 26224524). Sheep reactivity reported in scientific literature (PMID: 26563105).
Immunogen	A portion of amino acids 305-338 of human LOX protein were used as the immunogen.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation, Knockout Validated
Recommended Dilutions	Western Blot 1:100-1:2000, Immunohistochemistry reported in scientific literature (PMID 29404848), Immunocytochemistry/ Immunofluorescence 1:10-1:500, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Paraffin 1:10-1:500, Immunohistochemistry-Frozen, Knockout Validated reported in scientific literature (PMID 31173582)




Images


Western Blot: LOX Antibody [NBP2-24877] - Analysis of mouse ear lysate when LOX antibody (NBP2-24877) was used at 1:1000. The band shown is the expected inactive 50 kD proenzyme.




Immunocytochemistry/Immunofluorescence: LOX Antibody [NBP2-24877] - HeLa cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X TBS + 0.5% Triton-X100. The cells were incubated with anti-LOX (NBP2-24877) at a 1:200 dilution overnight at 4C and detected with an anti-rabbit Dylight 488 (Green) at a 1:500 dilution. Alpha tubulin was used as a co-stain at a 1:1000 dilution and detected with an anti-mouse Dylight 550 (Red) at a 1:500 dilution. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 40X objective.



Immunohistochemistry-Paraffin: LOX Antibody [NBP2-24877] - Analysis of LOX in formalin-fixed, paraffin-embedded mouse skin using an isotype control (top) and NBP2-24877 (bottom) at 1:100.



Western Blot: LOX Antibody [NBP2-24877] - Analysis of LOX in A) human skin and B) mouse ear skin lysate using NBP2-24877 at 1:500. Goat anti-rabbit IgG HRP secondary antibody and PicoTest ECL substrate solution were used for this test.



Publications

Ma, Y;Harris, J;Li, P;Jiang, C;Sun, H;Cao, H; An Integrative Transcriptome Subtraction Strategy to Identify Human lncRNAs That Specifically Play a Role in Activation of Human Hepatic Stellate Cells Non-coding RNA 2024-06-06 [PMID: 38921831]

Fuqua JD, Lawrence MM, Hettinger ZR et al. Impaired proteostatic mechanisms other than decreased protein synthesis limit old skeletal muscle recovery after disuse atrophy Journal of cachexia, sarcopenia and muscle 2023-07-14 [PMID: 37448295] (WB, Rat)

Wang C Regulatory Roles of Fibril-Forming Collagens in Cartilaginous Tissue Biomechanics and Mechanobiology Thesis 2023-01-01

Abbott C, Lawrence M, Kobak K et al. A novel stable isotope approach demonstrates surprising degree of age-related decline in skeletal muscle collagen proteostasis Function 2021-05-10 [PMID: 34124684] (WB, Mouse)

Cellular and Extracellular Matrix Basis for Heterogeneity in Mitral Annular Contraction Stephens EH, Fahrenholtz MM, Connell PS Cardiovasc Eng Technol [PMID: 26195991] (IHC-P, Sheep)

Wang C, Brisson BK, Terajima M et al. Type III Collagen is a Key Regulator of the Collagen Fibrillar Structure and Biomechanics of Articular Cartilage and Meniscus Matrix Biol. 2019-10-23 [PMID: 31655293] (WB, Mouse)

Terajima M, Taga Y, Cabral WA et al. Cyclophilin B control of lysine post-translational modifications of skin type I collagen PLoS Genet. 2019-06-01 [PMID: 31173582] (WB, IF/IHC, Knockout Validated, Mouse)

Rossow L, Veitl S, Vorlova S et al LOX-catalyzed collagen stabilization is a proximal cause for intrinsic resistance to chemotherapy. Oncogene. 2018-09-01 [PMID: 29780168] (Human)

Details:

Citation using the HRP version of this antibody.

Connell PS, Vekilov DP, Diaz CM et al. Eliminating Regurgitation Reduces Fibrotic Remodeling of Functional Mitral Regurgitation Conditioned Valves. Ann Biomed Eng 2018-02-05 [PMID: 29404848] (IF/IHC, Porcine)

Smith LR, Hammers DW, Sweeney HL, Barton ER. Increased collagen cross-linking is a signature of dystrophin-deficient muscle. Muscle Nerve. 2015-11-30 [PMID: 26616495] (WB, Mouse)

Stephens EH, Nguyen TC, Blazejewski JG et al. Extracellular matrix remodeling in wound healing of critical size defects in the mitral valve leaflet. Heart Vessels. 2015-11-12 [PMID: 26563105] (IHC-P, Sheep)

Connell PS, Azimuddin AF, Kim SE et al. Regurgitation Hemodynamics Alone Cause Mitral Valve Remodeling Characteristic of Clinical Disease States In Vitro. Ann Biomed Eng. 2015-07-30 [PMID: 26224524] (IHC-P, Porcine)

Details:

LOX antibody (Imgenex: IMG-6442A) was used for IHC-P analysis of lysyl oxidase protein in experiments involving organ cultures of mitral valves that were obtained from adolescent (6–9 months) porcine animals.

More publications at <http://www.novusbio.com/NBP2-24877>





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-24877

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB100-2527PEP	LOX Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-24877

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

