

Product Datasheet

SOD1/Cu-Zn SOD Antibody - BSA Free NBP2-24915

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 3 Publications: 21

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-24915

Updated 2/23/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-24915



NBP2-24915

SOD1/Cu-Zn SOD Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS


Product Description	
Host	Rabbit
Gene ID	6647
Gene Symbol	SOD1
Species	Human, Mouse, Rat, Bovine
Immunogen	The superoxide dismutase enzyme from bovine erythrocytes (CAS # 9054-89-1) was used as the immunogen.

Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1 - 2 ug/ml, Immunohistochemistry 1:200, Immunohistochemistry-Paraffin 1:200




Images


Western Blot: SOD1/Cu-Zn SOD Antibody [NBP2-24915] - Bovine ocular fluid. Antibody at 1:1000. Western blot image submitted by a verified customer review.




Immunohistochemistry-Paraffin: SOD1/Cu-Zn SOD Antibody [NBP2-24915] - Analysis of Superoxide Dismutase 1 in FFPE human liver tissue using an isotype control (top) and NBP2-24915 (bottom) at 5 ug/ml.



Western Blot: SOD1/Cu-Zn SOD Antibody [NBP2-24915] - Analysis of SOD in liver lysate of 1) human, 2) mouse, and 3) rat using this antibody. 0.25 ug/ml.



Western Blot: SOD1/Cu-Zn SOD Antibody [NBP2-24915] - Human ocular fluid. Antibody at 1:1000. Western blot image submitted by a verified customer review.



Publications

Mueller, J;van Muilekom, DR;Ehlers, J;Suhr, M;Hornburg, SC;Bang, C;Wilkes, M;Schultheiß, T;Maser, E;Rebl, A;Goldammer, T;Seibel, H;Schulz, C; Dietary Chlorella vulgaris supplementation modulates health, microbiota and the response to oxidative stress of Atlantic salmon Scientific reports 2024-10-10 [PMID: 39389986]

Bertolo A, Capossela S, Frankl G et al. Oxidative status predicts quality in human mesenchymal stem cells Stem Cell Res Ther Jan 6 2017 12:00AM [PMID: 28061861] (FLOW, Human). Stem Cell Res Ther. 2017-01-06 [PMID: 28061861]

Details:

This citation used the FITC version of this antibody.

Kwon I, Song W, Jang Y et al. Elevation of hepatic autophagy and antioxidative capacity by endurance exercise is associated with suppression of apoptosis in mice Annals of Hepatology [PMID: 31611063]

Chang YJ, Jenny LA, Li YS et al. CRISPR editing demonstrates rs10490924 raised oxidative stress in iPSC-derived retinal cells from patients with ARMS2/HTRA1-related AMD Proceedings of the National Academy of Sciences of the United States of America 2023-05-09 [PMID: 37126685]

F Wei, CJ Neal, TS Sakthivel, Y Fu, M Omer, A Adhikary, S Ward, KM Ta, S Moxon, M Molinari, J Asiatico, M Kinzel, SN Yarmolenko, V San Cheong, N Orlovskaya, R Ghosh, S Seal, M Coathup A novel approach for the prevention of ionizing radiation-induced bone loss using a designer multifunctional cerium oxide nanozyme Bioactive materials, 2022-09-21;21(0):547-565. 2022-09-21 [PMID: 36185749]

Noor MI, Rahman MS Roundup® disrupts tissue architecture, attenuates Na⁺/K⁺-ATPase expression, and induces protein oxidation/nitration, cellular apoptosis, and antioxidant enzyme expressions in the gills of goldfish, Carassius auratus Comparative biochemistry and physiology. Toxicology & pharmacology : CBP 2023-07-31 [PMID: 37532112]

Najjar R Raspberry Polyphenols Target Molecular Pathways of Heart Failure Thesis 2023-01-01

Martinez-Martel I, Bai X, Batalle G, Pol O New Treatment for the Cognitive and Emotional Deficits Linked with Paclitaxel-Induced Peripheral Neuropathy in Mice Antioxidants (Basel, Switzerland) 2022-12-01 [PMID: 36552595] (WB, Mouse)

Lacy B, Rahman MS Interactive effects of high temperature and pesticide exposure on oxidative status, apoptosis, and renin expression in kidney of goldfish: Molecular and cellular mechanisms of widespread kidney damage and renin attenuation Journal of applied toxicology : JAT 2022-06-13 [PMID: 35698815]

Lacy B, Rahman MS, Rahman MS Potential mechanisms of Na⁺/K⁺-ATPase attenuation by heat and pesticides co-exposure in goldfish: role of cellular apoptosis, oxidative/nitrative stress, and antioxidants in gills Environmental science and pollution research international [PMID: 35352221]

Carmona-Rodriguez L, Martinez-Rey D, Martin-Gonzalez P et al. Superoxide Dismutase-3 Downregulates Laminin alpha5 Expression in Tumor Endothelial Cells via the Inhibition of Nuclear Factor Kappa B Signaling Cancers 2022-02-26 [PMID: 35267534] (WB, Human)

Pouso-Vazquez E, Bai X, Batalle G et al. Effects of heme oxygenase 1 in the molecular changes and neuropathy associated with type 2 diabetes in mice Biochemical pharmacology 2022-03-08 [PMID: 35276215] (WB, Mouse)

More publications at <http://www.novusbio.com/NBP2-24915>



Procedures

Western Blot Protocol for SOD1/Cu-Zn SOD Antibody (NBP2-24915)

Western Blot Protocol

1. Perform SDS-PAGE on samples to be analyzed, loading 10-25 ug of total protein per lane.
2. Transfer proteins to PVDF membrane according to the instructions provided by the manufacturer of the membrane and transfer apparatus.
3. Stain the membrane with Ponceau S (or similar product) to assess transfer success, and mark molecular weight standards where appropriate.
4. Rinse the blot TBS -0.05% Tween 20 (TBST).
5. Block the membrane in 5% Non-fat milk in TBST (blocking buffer) for at least 1 hour.
6. Wash the membrane in TBST three times for 10 minutes each.
7. Dilute primary antibody in blocking buffer and incubate overnight at 4C with gentle rocking.
8. Wash the membrane in TBST three times for 10 minutes each.
9. Incubate the membrane in diluted HRP conjugated secondary antibody in blocking buffer (as per manufacturer's instructions) for 1 hour at room temperature.
10. Wash the blot in TBST three times for 10 minutes each (this step can be repeated as required to reduce background).
11. Apply the detection reagent of choice in accordance with the manufacturer's instructions.

Immunohistochemistry-Paraffin Protocol for SOD1/Cu-Zn SOD Antibody (NBP2-24915)

Immunohistochemistry-Paraffin Embedded Sections

Antigen Unmasking:

Bring slides to a boil in 10 mM sodium citrate buffer (pH 6.0) then maintain at a sub-boiling temperature for 10 minutes. Cool slides on bench-top for 30 minutes (keep slides in the sodium citrate buffer at all times).

Staining:

1. Wash sections in deionized water three times for 5 minutes each.
2. Wash sections in PBS for 5 minutes.
3. Block each section with 100-400 ul blocking solution (1% BSA in PBS) for 1 hour at room temperature.
4. Remove blocking solution and add 100-400 ul diluted primary antibody. Incubate overnight at 4 C.
5. Remove antibody solution and wash sections in wash buffer three times for 5 minutes each.
6. Add 100-400 ul HRP polymer conjugated secondary antibody. Incubate 30 minutes at room temperature.
7. Wash sections three times in wash buffer for 5 minutes each.
8. Add 100-400 ul DAB substrate to each section and monitor staining closely.
9. As soon as the sections develop, immerse slides in deionized water.
10. Counterstain sections in hematoxylin.
11. Wash sections in deionized water two times for 5 minutes each.
12. Dehydrate sections.
13. Mount coverslips.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-24915

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-34940-100ug	Recombinant Human SOD1/Cu-Zn SOD Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-24915

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

