

# Product Datasheet

## TRAILR1/TNFRSF10A Antibody [DyLight 755] NBP2-25092IR

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-25092IR](http://www.novusbio.com/NBP2-25092IR)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-25092IR](http://www.novusbio.com/reviews/destination/NBP2-25092IR)



**NBP2-25092IR**

TRAILR1/TNFRSF10A Antibody [DyLight 755]

| Product Information         |                                                                                                                                                                                                                                |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                         |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                    |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                       |
| Clonality                   | Polyclonal                                                                                                                                                                                                                     |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                             |
| Isotype                     | IgG                                                                                                                                                                                                                            |
| Conjugate                   | DyLight 755                                                                                                                                                                                                                    |
| Purity                      | Peptide affinity purified                                                                                                                                                                                                      |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                             |
| Product Description         |                                                                                                                                                                                                                                |
| Host                        | Rabbit                                                                                                                                                                                                                         |
| Gene ID                     | 8797                                                                                                                                                                                                                           |
| Gene Symbol                 | TNFRSF10A                                                                                                                                                                                                                      |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                              |
| Reactivity Notes            | 0                                                                                                                                                                                                                              |
| Specificity/Sensitivity     | Human TRAILR1/TNFRSF10A has only 1 isoform (468aa, 50kD). TRAILR1/TNFRSF10A antibody has no cross reaction to TRAILR2/TNFRSF10B.                                                                                               |
| Immunogen                   | Antibody was raised against a peptide corresponding to 19 amino acids near the carboxy terminus of human DR4 protein. The immunogen is located within the last 50 amino acids of DR4. Amino Acid Sequence: DSGKFIYLEDGTGSAVSLE |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                              |
| Product Application Details |                                                                                                                                                                                                                                |
| Applications                | Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                               |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                               |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                         |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP2-25092IR**

---

|               |                                               |
|---------------|-----------------------------------------------|
| NBP2-24891IR  | Rabbit IgG Isotype Control [DyLight 755]      |
| NBP2-33716PEP | TRAILR1/TNFRSF10A Recombinant Protein Antigen |
| 210-TA-005    | TNF-alpha [Unconjugated]                      |
| 6949-DR-100   | TRAILR1/TNFRSF10A [Unconjugated]              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-25092IR](http://www.novusbio.com/reviews/submit/NBP2-25092IR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

