

Product Datasheet

Inhibin alpha Antibody NBP2-26096

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-26096

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-26096



NBP2-26096**Inhibin alpha Antibody****Product Information**

Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA

Product Description


Description	Novus Biologicals Goat Inhibin alpha Antibody (NBP2-26096) is a polyclonal antibody validated for use in IHC, WB and ELISA. Anti-Inhibin alpha Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Goat
Gene ID	3623
Gene Symbol	INHA
Species	Human, Bovine
Immunogen	Peptide with sequence C-RPLHVRTTSDGGYS corresponding to C-Terminus according to NP_002182.1.

Product Application Details

Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Peptide ELISA
Recommended Dilutions	Western Blot 1:100 - 1:2000, Immunohistochemistry, Immunohistochemistry-Paraffin, Peptide ELISA Detection limit 1:128000
Application Notes	This Inhibin alpha antibody is validated for IHC-P from a verified customer review.

Images

Immunohistochemistry-Paraffin: Inhibin alpha Antibody [NBP2-26096] - 60X image of Hu Testis stained with Inhibin alpha antibody at 1:50 shows questionable punctate or vesicular staining, both cytoplasmic and somewhat nuclear, in fluorescent green 488. Did not correspond completely with multiple other INHA antibodies tested. IHC staining with citrate pH 6 antigen retrieval; detected staining with a 488-conjugated anti-goat secondary antibody. IHC image submitted by a verified customer review.

**Publications**

Nistal M, Pastrian LG, Gonzalez-Peramato P, De Miguel MP. Inhibin bodies: a new marker for immature Sertoli cells. Histopathology 2011-06-01 [PMID: 21707703]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-26096

HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control
NBP1-87563PEP	Inhibin alpha Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-26096

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



