

Product Datasheet

IL-23R Antibody (15N6C6) - BSA Free NBP2-27091

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-27091

Updated 2/23/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-27091



NBP2-27091

IL-23R Antibody (15N6C6) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	15N6C6
Preservative	0.05% Sodium Azide
Isotype	IgG1 Kappa
Purity	Protein G purified
Buffer	PBS
Product Description	
Host	Mouse
Gene ID	149233
Gene Symbol	IL23R
Species	Human, Primate
Reactivity Notes	Chimpanzee, Orangutan (100%), Mouse, Porcine, Canine, Rhesus monkey (89%), Rat, Equine, Bovine (84%)
Immunogen	A synthetic peptide corresponding to amino acids 288-315 of hIL-23 receptor was used as immunogen for this antibody.
Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1 - 3 ug/mL, Immunohistochemistry, Immunohistochemistry-Paraffin
Application Notes	There are multiple isoforms of IL-23 receptor. Researchers are encouraged to use a tool like BLAST to determine the suitability of this antibody for their testing.




Images

Western Blot: IL-23 R Antibody (15N6C6) [NBP2-27091] - Analysis using IL-23 receptor antibody. Lysate from (A) human HCT116, (B) Jurkat and C) mouse NIH 3T3 cells probed with IL-23 receptor antibody at 1 ug/mL.



Immunohistochemistry-Paraffin: IL-23 R Antibody (15N6C6) [NBP2-27091] - Human Appendix FFPE tissue section stained with IL-23 R antibody. Image from verified customer review.



Immunohistochemistry-Paraffin: IL-23 R Antibody (15N6C6) [NBP2-27091] - Analysis of IL-23 receptor in human colon tissue using IL-23 receptor antibody at 5 ug/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-27091

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
NBP2-27092B	IL-23R Antibody (15N6C6) [Biotin]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-27091

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

