

Product Datasheet

Dnmt3L Antibody NBP2-27099

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-27099

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-27099



NBP2-27099


Dnmt3L Antibody

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	0.2 ml PBS and 0.05% BSA
Product Description	
Host	Rabbit
Gene ID	29947
Gene Symbol	DNMT3L
Species	Human, Primate
Reactivity Notes	The amino acid sequence used as immunogen for this antibody is 100% homologous in White-cheeked Gibbon, 95% in Macaque, 89% in Bovine, New World Monkey, and Equine, and 84% in Porcine, and Opossum.
Immunogen	A portion of amino acids 25-75 from human Dnmt3L was used as the immunogen for this antibody.
Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 2-4ug/ml~, Immunohistochemistry, Immunohistochemistry-Paraffin 5ug/ml




Images

Western Blot: Dnmt3L Antibody [NBP2-27099] - Analysis of Dnmt3L using Dnmt3L antibody. GST tagged full length recombinant Human Dnmt3L in the 1) pre immune sera and 2) absence of immunizing peptide, probed with 2 ug/ml of Dnmt3L antibody. Goat anti-rabbit Ig HRP secondary antibody and PicoTect ECL substrate solution were used for this test.



Immunohistochemistry-Paraffin: Dnmt3L Antibody [NBP2-27099] - Analysis of Dnmt3L in human kidney using Dnmt3L antibody at 5 ug/ml.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-27099

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-51647-0.02mg	Recombinant Human Dnmt3L His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-27099

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

