

Product Datasheet

IL-17/IL-17A Antibody (4H1524) NBP2-27338

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-27338

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-27338



NBP2-27338


IL-17/IL-17A Antibody (4H1524)

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4H1524
Preservative	0.05% Sodium Azide
Isotype	IgG2b Kappa
Purity	Protein G purified
Buffer	0.2 ml PBS and 0.05% BSA
Product Description	
Host	Mouse
Gene ID	3605
Gene Symbol	IL17A
Species	Human
Immunogen	A portion of amino acids 1-100 was used as immunogen
Product Application Details	
Applications	Western Blot, Flow Cytometry, Block/Neutralize
Recommended Dilutions	Western Blot 1-3ug/ml~, Flow Cytometry 1:10-1:1000, Block/Neutralize




Images


Western Blot: IL-17/IL-17A Antibody (4H1524) [NBP2-27338] - Analysis of IL-17A in 1) human, 2) mouse and 3) rat full-length recombinant IL-17A protein using NBP2-27338 at 1 ug/ml.



Flow Cytometry: IL-17/IL-17A Antibody (4H1524) [NBP2-27338] - Analysis using the PE conjugate of NBP2-25229. IL-17A positive cells in PMA/IM stimulated fresh PBMC: PBMC stimulated for 6h with PMA (25ng/ml) & IM (250ng/ml), in the presence of BFA. Cells were surface stained with CD4 APC, fixed & permeabilised for intracellular staining with PE conjugated



Block/Neutralize: IL-17/IL-17A Antibody (4H1524) [NBP2-27338] - Analysis using Azide Free version of NBP2-27338. Neutralization by Human IL-17A Antibody. Secretion of IL-6 by 3T3 cells stimulated with recombinant human IL-17A (20ng/ml) is neutralized by increasing concentrations of IL-17A monoclonal antibody. Antibody



Publications

Liu G, Wu H, Chen L et al. Effects of interleukin-17 on human retinal vascular endothelial cell capillary tube formation in vitro. Mol Med Rep. 2017-07-01 [PMID: 28560397]

Wen W, Wan Z, Ren K et al. Potassium supplementation inhibits IL-17A production induced by salt loading in human T lymphocytes via p38/MAPK-SGK1 pathway. Experimental and Molecular Pathology. 2016-03-25 [PMID: 27020669] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-27338

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
NBP2-27252	IL-17/IL-17A Antibody (4H1524) [PE]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-27338

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

