

Product Datasheet

ER beta/NR3A2 Antibody (ERb455) NBP2-29439

Unit Size: 0.1 mg

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-29439

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-29439



NBP2-29439

ER beta/NR3A2 Antibody (ERb455)

Product Information	
Unit Size	0.1 mg
Concentration	0.2 mg/ml
Storage	Store at 4C.
Clonality	Monoclonal
Clone	ERb455
Preservative	0.05% Sodium Azide
Isotype	IgG2a Kappa
Purity	Protein A or G purified
Buffer	10 mM PBS with 0.05% BSA


Product Description	
Description	200ug/ml of antibody purified from Bioreactor Concentrate by Protein A or G. Prepared in 10 mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0 mg/ml. (NBP2-33066) Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80C.
Host	Mouse
Gene ID	2100
Gene Symbol	ESR2
Species	Human, Mouse, Rat, Porcine, Equine, Monkey, Sheep
Immunogen	C-terminus fragment of recombinant human ER beta/NR3A2 protein (Uniprot: Q92731)

Product Application Details	
Applications	Western Blot, Simple Western, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1-2 ug/ml, Simple Western, Flow Cytometry 1-2 ug/million cells, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1-2 ug/ml, Immunohistochemistry-Paraffin 1-2 ug/ml
Application Notes	Immunohistochemistry (Formalin-fixed): 1-2ug/ml for 30 minutes at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes. Optimal dilution for a specific application should be determined.




Images


Western Blot: ER beta/NR3A2 Antibody (ERb455) [NBP2-29439] - Western Blot Analysis of human MCF-7 cell lysate using ER beta/NR3A2 Antibody (ERb455).




Immunohistochemistry-Paraffin: ER beta/NR3A2 Antibody (ERb455) [NBP2-29439] - Formalin-fixed, paraffin-embedded human Bladder Carcinoma stained with ER beta/NR3A2 Antibody (ERb455).




Flow Cytometry: ER beta/NR3A2 Antibody (ERb455) [NBP2-29439] - Flow Cytometry of human ER beta on BT474 cells. Black: cells alone; Grey: Isotype Control; Green: AF488-labeled ER beta/NR3A2 Mouse Monoclonal Antibody (ERb455).




Immunohistochemistry-Paraffin: ER beta/NR3A2 Antibody (ERb455) [NBP2-29439] - Formalin paraffin breast cancer stained with Estrogen receptor beta Ab (ERb455).



Immunohistochemistry-Paraffin: ER beta/NR3A2 Antibody (ERb455) [NBP2-29439] - Formalin-paraffin ovarian cancer stained with Estrogen receptor beta Ab (ERb455).




Immunohistochemistry-Paraffin: ER beta/NR3A2 Antibody (ERb455) [NBP2-29439] - Formalin-fixed, paraffin-embedded human Ovarian Carcinoma stained with ER beta/NR3A2 Antibody (ERb455).



Immunohistochemistry-Paraffin: ER beta/NR3A2 Antibody (ERb455) [NBP2-29439] - Formalin-fixed, paraffin-embedded human Gastric Carcinoma stained with ER beta/NR3A2 Antibody (ERb455).



Flow Cytometry: ER beta/NR3A2 Antibody (ERb455) [NBP2-29439] - Flow Cytometry for human ER-beta on MCF-7 cells. Black: cells alone; Green: Isotype Control; Red: PE-labeled ER beta/NR3A2 Antibody (ERb455).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-29439

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
H00002100-Q01-10ug	Recombinant Human ER beta/NR3A2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-29439

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

