

Product Datasheet

Olig2 [p Ser14, p Ser10, p Ser13] Antibody - Azide Free NBP2-29523

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-29523

Updated 2/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-29523




NBP2-29523

Olig2 [p Ser14, p Ser10, p Ser13] Antibody - Azide Free

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	10mM HEPES (pH 7.5), 0.15M NaCl, 0.1 mg/ml BSA and 50% Glycerol
Target Molecular Weight	32 kDa
Product Description	
Host	Rabbit
Gene ID	10215
Gene Symbol	OLIG2
Species	Mouse, Rat
Marker	Oligodendrocyte Marker
Specificity/Sensitivity	Specific for endogenous levels of the ~32 kDa Olig2 phosphorylated at Ser10, Ser 13, and Ser14. Immunolabeling is blocked by preadsorption with the phosphopeptide used as antigen, but not by the corresponding non-phosphopeptide.
Immunogen	Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser10,13,14 conjugated to KLH. Accession # Q13516
Preparation Method	Prepared from rabbit serum by affinity purification via sequential chromatography on phospho- and dephosphopeptide affinity columns.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:1000

Images

Western Blot: OLIG2 [p Ser14, p Ser10, p Ser13] Antibody [NBP2-29523] - Western blot of rat neonatal brain lysate showing specific immunolabeling of the ~ 32k form of the Olig2 protein phosphorylated at Ser10,13,14. Immunolabeling is blocked by the phospho-peptide used as antigen (peptide), but not by the corresponding dephospho-peptide (not shown).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-29523

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00010215-P01-10ug	Recombinant Human Olig2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-29523

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

