

# Product Datasheet

## UBE4A Antibody [FITC] NBP2-32200F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-32200F](http://www.novusbio.com/NBP2-32200F)

Updated 11/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-32200F](http://www.novusbio.com/reviews/destination/NBP2-32200F)



**NBP2-32200F**

UBE4A Antibody [FITC]

| Product Information         |                                                                                                                                                                                            |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                     |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                   |
| Clonality                   | Polyclonal                                                                                                                                                                                 |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                         |
| Isotype                     | IgG                                                                                                                                                                                        |
| Conjugate                   | FITC                                                                                                                                                                                       |
| Purity                      | Immunogen affinity purified                                                                                                                                                                |
| Buffer                      | PBS                                                                                                                                                                                        |
| Product Description         |                                                                                                                                                                                            |
| Host                        | Rabbit                                                                                                                                                                                     |
| Gene ID                     | 9354                                                                                                                                                                                       |
| Gene Symbol                 | UBE4A                                                                                                                                                                                      |
| Species                     | Human                                                                                                                                                                                      |
| Immunogen                   | The immunogen this antibody was made to, maps to a region between residue 1016 to 1066 of human Ubiquitination Factor E4A using the numbering given in entry NP_001191006.1 (GeneID 9354). |
| Product Application Details |                                                                                                                                                                                            |
| Applications                | Western Blot, Immunoprecipitation                                                                                                                                                          |
| Recommended Dilutions       | Western Blot, Immunoprecipitation                                                                                                                                                          |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                     |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-32200F**

---

|                |                                           |
|----------------|-------------------------------------------|
| NBP2-24892     | Rabbit IgG Isotype Control [FITC]         |
| NBP2-47590PEP  | UBE4A Recombinant Protein Antigen         |
| 211-TBB-010/CF | Lymphotoxin-alpha/TNF-beta [Unconjugated] |
| NBP2-56653PEP  | UBE4A Recombinant Protein Antigen         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-32200F](http://www.novusbio.com/reviews/submit/NBP2-32200F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

