

# Product Datasheet

## Kinesin 5A Antibody [PE] NBP2-32231PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-32231PE](http://www.novusbio.com/NBP2-32231PE)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-32231PE](http://www.novusbio.com/reviews/destination/NBP2-32231PE)



**NBP2-32231PE**

Kinesin 5A Antibody [PE]

| Product Information         |                                                                                                                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                           |
| Storage                     | Store at 4C in the dark.                                                                                                                                                              |
| Clonality                   | Polyclonal                                                                                                                                                                            |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                    |
| Isotype                     | IgG                                                                                                                                                                                   |
| Conjugate                   | PE                                                                                                                                                                                    |
| Purity                      | Immunogen affinity purified                                                                                                                                                           |
| Buffer                      | PBS                                                                                                                                                                                   |
| Product Description         |                                                                                                                                                                                       |
| Host                        | Rabbit                                                                                                                                                                                |
| Gene ID                     | 3798                                                                                                                                                                                  |
| Gene Symbol                 | KIF5A                                                                                                                                                                                 |
| Species                     | Human                                                                                                                                                                                 |
| Immunogen                   | The immunogen this antibody was made to, maps to a region between residue 982 to 1032 of human Kinesin Family Member 5A using the numbering given in entry NP_004975.2 (GeneID 3798). |
| Product Application Details |                                                                                                                                                                                       |
| Applications                | Western Blot, Immunoprecipitation                                                                                                                                                     |
| Recommended Dilutions       | Western Blot, Immunoprecipitation                                                                                                                                                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-32231PE**

---

|               |                                           |
|---------------|-------------------------------------------|
| NBP2-24983    | Rabbit IgG Isotype Control [PE]           |
| NBP2-57903PEP | Kinesin 5A Recombinant Protein Antigen    |
| NBP3-39712    | Human Kinesin 5A ELISA Kit (Colorimetric) |
| NB500-181     | Eg5 Antibody - Azide and BSA Free         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-32231PE](http://www.novusbio.com/reviews/submit/NBP2-32231PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

