

Product Datasheet

CD45RO Antibody (UCHL-1) - Azide and BSA Free NBP2-33104

Unit Size: 0.1 mg

Store at -20 to -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-33104

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-33104



NBP2-33104

CD45RO Antibody (UCHL-1) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at -20 to -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	UCHL-1
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Protein A or G purified
Buffer	10 mM PBS
Target Molecular Weight	182 kDa


Product Description	
Description	1.0 mg/ml of antibody purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS WITHOUT BSA & azide. Also available at 200 ug/ml WITH BSA & azide (NBP2-29631). Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80C.
Host	Mouse
Gene ID	5788
Gene Symbol	PTPRC
Species	Human, Chimpanzee, Monkey
Reactivity Notes	Common Marmoset
Marker	T-Cell Marker
Specificity/Sensitivity	Recognizes a 180-185kDa protein, identified as isoform of leukocyte common antigen (CD45RO) (4th Leucocyte Typing Workshop: Code No. N31). The epitope recognized by this antibody is sensitive to neuraminidase digestion. This antibody reacts with mature activated T-cells, most thymocytes, and a sub-population of resting T-cells within both CD4 and CD8 subsets. It shows no reactivity with normal B or natural killer cells, but reacts with granulocytes and monocytes. Reportedly, it is useful to identify T-cell lymphomas and leukemias. It rarely stains NK cells or B-cell lymphomas.
Immunogen	Cultured human T-cells from an IL-2-dependent T-cell line (CA1)

Product Application Details	
Applications	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready, Dual RNAscope ISH-IHC
Recommended Dilutions	Flow Cytometry 0.5-1ug/million cells, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 0.5-1ug/ml, Immunohistochemistry-Paraffin 0.5-1.0ug/ml, CyTOF-ready, Dual RNAscope ISH-IHC
Application Notes	Immunohistochemistry (Formalin-fixed): 0.25-0.5ug/ml for 30 minutes at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes. Optimal dilution for a specific application should be determined.




Images


Immunocytochemistry/Immunofluorescence: CD45RO Antibody (UCHL-1) - Azide and BSA Free [NBP2-33104] - Immunofluorescence Analysis of human tonsil FFPE section stained with CD45RO antibody (UCHL-1) and goat anti-mouse HRP, detected using CF488A tyramide in Tyramide Amplification Buffer Plus.




Immunohistochemistry-Paraffin: CD45RO Antibody (UCHL-1) - Azide and BSA Free [NBP2-33104] - Staining of human tonsil




Flow Cytometry: CD45RO Antibody (UCHL-1) - Azide and BSA Free [NBP2-33104] - CD45RO expression by lymphocyte gated population of PBMC: PBMC were stained either with purified CD45RO antibody (0.5 ug/10⁶ cells) or isotype control following surface staining protocol. PE-conjugated goat anti-mouse IgG, secondary antibody used to identify the positive staining on PPI negative lymphocyte gated population.



Flow Cytometry: CD45RO Antibody (UCHL-1) - Azide and BSA Free [NBP2-33104] - Flow Cytometric Analysis of PFA-fixed Jurkat cells using CD45RO Antibody (UCHL-1) followed by Goat anti- Mouse- IgG-CF488 (Blue); Isotype Control (Red).



Dual RNAscope ISH-IHC: CD45RO Antibody (UCHL-1) - Azide and BSA Free [NBP2-33104] - Formalin-fixed paraffin-embedded tissue sections of human tonsil were probed for CD45 mRNA (ACD RNAScope probe, catalog # 601991; Fast Red chromogen, ACD catalog # 322500) followed by immunohistochemistry using mouse monoclonal (Novus catalog # NBP2-29631) at 5ug/mL for 1 hour at room temperature and detection with anti-mouse IgG VisUCyte HRP Polymer Antibody (R&D Systems Catalog # VC001) and green HRP chromogen (ACD Catalog # 322500).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-33104

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
NBP2-33104AF405	CD45RO Antibody (UCHL-1) [Alexa Fluor® 405]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-33104

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

