

# Product Datasheet

## Maxi Potassium channel alpha Antibody - BSA Free NBP2-33726

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-33726](http://www.novusbio.com/NBP2-33726)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-33726](http://www.novusbio.com/reviews/destination/NBP2-33726)



**NBP2-33726**

Maxi Potassium channel alpha Antibody - BSA Free

**Product Information**

|                      |                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                                           |
| <b>Concentration</b> | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| <b>Clonality</b>     | Polyclonal                                                                                                       |
| <b>Preservative</b>  | 0.02% Sodium Azide                                                                                               |
| <b>Isotype</b>       | IgG                                                                                                              |
| <b>Purity</b>        | Immunogen affinity purified                                                                                      |
| <b>Buffer</b>        | PBS (pH 7.2) and 40% Glycerol                                                                                    |

**Product Description**


|                    |                                                                                                                                                                                                                     |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Novus Biologicals Rabbit Maxi Potassium channel alpha Antibody - BSA Free (NBP2-33726) is a polyclonal antibody validated for use in IHC. All Novus Biologicals antibodies are covered by our 100% guarantee.       |
| <b>Host</b>        | Rabbit                                                                                                                                                                                                              |
| <b>Gene ID</b>     | 3778                                                                                                                                                                                                                |
| <b>Gene Symbol</b> | KCNMA1                                                                                                                                                                                                              |
| <b>Species</b>     | Human                                                                                                                                                                                                               |
| <b>Immunogen</b>   | This antibody was developed against a recombinant protein corresponding to amino acids:<br>INLCDMCVILSANQNNIDDTSLQDKECILASLNIKSMQFDDSIGVLQANSQGFTP<br>PGMDRSPDNPVHGMLRQPSITTGVNIPITELVNDTNVQFLDQDDDDDPDTE<br>LYLTQP |

**Product Application Details**


|                              |                                                                                 |
|------------------------------|---------------------------------------------------------------------------------|
| <b>Applications</b>          | Immunohistochemistry-Paraffin, Immunohistochemistry                             |
| <b>Recommended Dilutions</b> | Immunohistochemistry 1:200 - 1:500, Immunohistochemistry-Paraffin 1:200 - 1:500 |
| <b>Application Notes</b>     | For IHC-Paraffin, HIER pH 6 retrieval is recommended.                           |

**Images**


Immunohistochemistry-Paraffin: Maxi Potassium channel alpha Antibody [NBP2-33726] - Staining of human Fallopian tube shows weak membranous positivity in glandular cells.




Immunohistochemistry-Paraffin: Maxi Potassium channel alpha Antibody [NBP2-33726] - Staining of human rectum shows moderate membranous positivity in glandular cells.




Immunohistochemistry-Paraffin: Maxi Potassium channel alpha Antibody [NBP2-33726] - Staining of human fallopian tube shows weak membranous positivity in glandular cells.



Immunohistochemistry-Paraffin: Maxi Potassium channel alpha Antibody [NBP2-33726] - Staining of human skeletal muscle shows no positivity in myocytes as expected.



Immunohistochemistry-Paraffin: Maxi Potassium channel alpha Antibody [NBP2-33726] - Staining of human cerebellum shows moderate membranous positivity in purkinje cells.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-33726**

---

|               |                                                          |
|---------------|----------------------------------------------------------|
| NBP2-33726PEP | Maxi Potassium channel alpha Recombinant Protein Antigen |
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]            |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]      |
| NBP2-24891    | Rabbit IgG Isotype Control                               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-33726](http://www.novusbio.com/reviews/submit/NBP2-33726)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

