

Product Datasheet

CD37 Antibody - BSA Free NBP2-33970

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-33970

Updated 2/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-33970



NBP2-33970

CD37 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Host	Rabbit
Gene ID	951
Gene Symbol	CD37
Species	Human
Immunogen	This antibody was developed against a recombinant protein corresponding to amino acids: ISTQRAQLERSLRDVVEKTIQKYGTNPEETAEEESWDYVQFQLRCCGWHYPQ DWFQVLILRGNGSEA


Product Application Details	
Applications	Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry 1:50 - 1:200, Immunohistochemistry-Paraffin 1:50 - 1:200
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining of human tonsil shows strong membranous positivity in non-germinal center cells.




Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining of human cerebral cortex shows low expression as expected.




Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining of human tonsil shows high expression.




Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining in human tonsil and cerebral cortex tissues using anti-CD37 antibody. Corresponding CD37 RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining of human lymph node.




Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining of human testis.



Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining of human spleen.




Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Analysis in human tonsil and skeletal muscle tissues. Corresponding CD37 RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining of human cerebral cortex shows no positivity in neurons as expected.



Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining of human cerebral cortex, lymph node, skeletal muscle and tonsil using Anti-CD37 antibody NBP2-33970 (A) shows similar protein distribution across tissues to independent antibody NBP2-33969 (B).



Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining of human lymph node shows strong membranous positivity in non-germinal center cells.



Immunohistochemistry-Paraffin: CD37 Antibody [NBP2-33970] - Staining of human skeletal muscle shows no positivity in myocytes as expected.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-33970

NBP2-33970PEP	CD37 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-33970

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

