

# Product Datasheet

## IL-37/IL-1F7 Antibody - BSA Free NBP2-34208

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-34208](http://www.novusbio.com/NBP2-34208)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-34208](http://www.novusbio.com/reviews/destination/NBP2-34208)



**NBP2-34208**

IL-37/IL-1F7 Antibody - BSA Free

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| Clonality           | Polyclonal                                                                                                       |
| Preservative        | 0.02% Sodium Azide                                                                                               |
| Isotype             | IgG                                                                                                              |
| Purity              | Immunogen affinity purified                                                                                      |
| Buffer              | PBS (pH 7.2) and 40% Glycerol                                                                                    |


| Product Description |                                                                                                                                                                                                          |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Novus Biologicals Rabbit IL-37/IL-1F7 Antibody - BSA Free (NBP2-34208) is a polyclonal antibody validated for use in IHC and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                | Rabbit                                                                                                                                                                                                   |
| Gene ID             | 27178                                                                                                                                                                                                    |
| Gene Symbol         | IL37                                                                                                                                                                                                     |
| Species             | Human                                                                                                                                                                                                    |
| Immunogen           | This antibody was developed against a recombinant protein corresponding to amino acids:<br>ENSGVKMGSEDWEKDEPQCCLEDPVSPLEPGPSLPAMNFBVHTSPKVKNLNP<br>KKFSIHDQDQHKVLVLDSGNLIAVPDKNYIRPEIFFALA               |

| Product Application Details |                                                                                                                                     |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry                                        |
| Recommended Dilutions       | Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200 |
| Application Notes           | For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.                       |




## Images

Immunocytochemistry/Immunofluorescence: IL-37/IL-1F7 Antibody [NBP2-34208] - Staining of human cell line PC-3 shows localization to nucleoplasm & vesicles. Antibody staining is shown in green.



Immunohistochemistry-Paraffin: IL-37/IL-1F7 Antibody [NBP2-34208] - Staining of human spleen shows positivity in white pulp.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-34208**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| NBP2-34208PEP | IL-37/IL-1F7 Recombinant Protein Antigen            |
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-34208](http://www.novusbio.com/reviews/submit/NBP2-34208)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

