

Product Datasheet

Mitochondria Antibody (113-1) - Azide and BSA Free NBP2-34517-0.1mg

Unit Size: 0.1 mg

Store at -20 to -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-34517

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-34517



NBP2-34517-0.1mg

Mitochondria Antibody (113-1) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at -20 to -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	113-1
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Protein A or G purified
Buffer	10 mM PBS
Target Molecular Weight	60 kDa


Product Description	
Description	1.0 mg/ml of antibody purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS WITHOUT BSA. Also available at 200 ug/ml WITH BSA & azide (NBP2-32980). Antibody with azide - store at 2 to 8C.
Host	Mouse
Species	Human, Mouse (Negative), Rat (Negative)
Reactivity Notes	Does not react with Mouse or Rat.
Marker	Marker for Human Cells, Granular RCC s & Salivary Tumor
Specificity/Sensitivity	This monoclonal antibody recognizes a 60kDa antigen associated with the mitochondria in human cells. It can be used to stain mitochondria in cell or tissue preparations and can be used as a mitochondrial marker in subcellular fractions. It produces a spaghetti-like pattern in normal and malignant cells. This antibody is an excellent marker for human cells in xenographic model research. It reacts specifically with human cells, including neurons and embryonic stem cells. Immunostaining pattern with anti-mitochondrial monoclonal antibody has been reported as a useful discriminatory adjunct in the complex differential diagnosis of granular renal cell tumors. Reportedly, this monoclonal antibody facilitates the classification of salivary tumors.
Immunogen	Semi-Purified mitochondrial preparation

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 0.25-0.5ug/ml, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 0.5-1ug/ml, Immunohistochemistry-Paraffin 0.5-1ug/ml, Immunohistochemistry-Frozen 0.5-1ug/mlimmunohistochemistry-Paraffin 0.5-1ug/ml
Application Notes	Immunohistochemistry (Formalin-fixed): 1-2ug/ml for 30 minutes at RT. Staining of formalin-fixed tissues is enhanced by heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes. Optimal dilution for a specific application should be determined.




Images


Western Blot: Mitochondria Antibody (113-1) - Azide and BSA Free [NBP2-34517] - Western Blot analysis of HeLa cell lysate using Mitochondria Antibody (113-1).




Immunocytochemistry/Immunofluorescence: Mitochondria Antibody (113-1) - Azide and BSA Free [NBP2-34517] - Confocal immunofluorescence image of HeLa cells using labeled is Green (CF488) and Red Dot is used to label the nuclei Red.




Immunohistochemistry-Paraffin: Mitochondria Antibody (113-1) - Azide and BSA Free [NBP2-34517] - Biotin conjugated Mitochondria expression in human prostate cancer tissue (left) and mouse stomach tissue (right). Image from verified customer review.




Immunohistochemistry-Paraffin: Mitochondria Antibody (113-1) - Azide and BSA Free [NBP2-34517] - Formalin-fixed paraffin-embedded human Renal cell carcinoma stained with Mitochondria Monoclonal antibody (113-1).




Immunohistochemistry-Paraffin: Mitochondria Antibody (113-1) - Azide and BSA Free [NBP2-34517] - Formalin-fixed, paraffin-embedded human testicular carcinoma stained with mitochondria monoclonal antibody (113-1).



Immunohistochemistry-Paraffin: Mitochondria Antibody (113-1) - Azide and BSA Free [NBP2-34517] - Formalin-fixed paraffin-embedded human pancreas stained with mitochondria monoclonal antibody (113-1).



Immunohistochemistry-Frozen: Mouse Monoclonal Mitochondria Antibody (113-1) - Azide and BSA Free [NBP2-34517] - H2030 cells (human lung tumor cell line) have been injected intracranially in immunodeficient mice pups to form tumors. NBP2-34517AF488 (green) detects human cells in murine brain. Nuclei have been stained with Hoechst (blue). Image from a verified customer review.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-34517-0.1mg

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
NBP2-29448	Mitochondria Isolation Kit

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-34517

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

