

# Product Datasheet

## **Chromogranin A Antibody (LK2H10 + PHE5 + CGA/414) [Janelia Fluor® 525] NBP2-34674JF525**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-34674JF525](http://www.novusbio.com/NBP2-34674JF525)

Updated 11/11/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-34674JF525](http://www.novusbio.com/reviews/destination/NBP2-34674JF525)



**NBP2-34674JF525**

Chromogranin A Antibody (LK2H10 + PHE5 + CGA/414) [Janelia Fluor® 525]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | LK2H10 + PHE5 + CGA/414                                                                     |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1 Kappa                                                                                  |
| Conjugate           | Janelia Fluor 525                                                                           |
| Purity              | Protein A or G purified                                                                     |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description |                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------|
| Host                | Mouse                                                                                                        |
| Gene ID             | 1113                                                                                                         |
| Gene Symbol         | CHGA                                                                                                         |
| Species             | Human, Mouse, Rat, Porcine, Monkey                                                                           |
| Marker              | Neuroendocrine Marker                                                                                        |
| Immunogen           | Human pheochromocytoma (LK2H10 & PHE5) & recombinant human chromogranin A protein (CGA414) (Uniprot: P10645) |
| Notes               | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.                        |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Immunohistochemistry-Paraffin, Immunohistochemistry                    |
| Recommended Dilutions       | Immunohistochemistry, Immunohistochemistry-Paraffin                    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-34674JF525**

---

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| NBP1-43319JF525 | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [Janelia Fluor 525] |
| NBC1-26369      | Recombinant Human Chromogranin A Protein                          |
| DKK300          | Kallikrein 3/PSA [HRP]                                            |
| DY9098-05       | Chromogranin A [Biotin]                                           |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-34674JF525](http://www.novusbio.com/reviews/submit/NBP2-34674JF525)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

