

Product Datasheet

Nucleoli Marker Antibody (NM95) - Azide and BSA Free NBP2-34695-0.1mg

Unit Size: 0.1 mg

Store at -20 to -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-34695

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-34695




NBP2-34695-0.1mg

Nucleoli Marker Antibody (NM95) - Azide and BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at -20 to -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	NM95
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Protein A or G purified
Buffer	10 mM PBS
Product Description	
Description	1.0 mg/ml of antibody purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS WITHOUT BSA & azide. Also available at 200 ug/ml WITH BSA & azide (NBP2-32886). Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80C.
Host	Mouse
Species	Human, Bovine (Negative), Mouse (Negative), Rat (Negative)
Reactivity Notes	Does not react with Mouse, Rat or Bovine.
Marker	Marker For Human Cells
Specificity/Sensitivity	This monoclonal antibody is an excellent marker for human cells in xenographic model research. It reacts specifically with human cells. This monoclonal antibody is part of a new panel of reagents, which recognizes subcellular organelles or compartments of human cells. These markers may be useful in identification of these organelles in cells, tissues, and biochemical preparations. monoclonal antibody NM95 recognizes an antigen associated with the nucleoli in human cells. It can be used to stain the nucleoli in cell or tissue preparations and can be used as a marker of the nucleoli in subcellular fractions. It produces a speckled pattern in the nuclei of cells of normal and malignant cells and may be used to stain the nucleoli of cells in fixed or frozen tissue sections. It can be used with paraformaldehyde fixed frozen tissue or cell preparations and formalin fixed, paraffin-embedded tissue sections.
Immunogen	Nuclei of myeloid leukemia biopsy cells
Product Application Details	
Applications	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin
Recommended Dilutions	Flow Cytometry, Immunocytochemistry/ Immunofluorescence 0.5-1.0ug/ml, Immunohistochemistry-Paraffin
Application Notes	Immunocytochemistry (Acetone-fixed cells): 1-2ug/ml for 30 minutes at RT. Immunohistochemistry (Formalin-fixed): 0.25-0.5ug/ml for 30 minutes at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes. Optimal dilution for a specific application should be determined.

Images


Immunocytochemistry/Immunofluorescence: Nucleoli Marker Antibody (NM95) - Azide and BSA Free [NBP2-34695] - Immunofluorescence Analysis of PFA-fixed K562 cells labeling Nucleolin with followed by Goat anti-Mouse IgG-CF488 (Green). Membrane is stained with Phalloidin-CF640.




Immunohistochemistry-Paraffin: Nucleoli Marker Antibody (NM95) - Azide and BSA Free [NBP2-34695] - Staining of human colon ca.



Flow Cytometry: Nucleoli Marker Antibody (NM95) - Azide and BSA Free [NBP2-34695] - Flow Cytometry of human nucleolar Ag on 293T cells
Black: cells alone: Grey: isotype control: Green: AF488-labeled human Nucleolar Monoclonal antibody (NM95)




Immunohistochemistry-Paraffin: Nucleoli Marker Antibody (NM95) - Azide and BSA Free [NBP2-34695] - Formalin-fixed, paraffin-embedded human Colon Carcinoma stained with Nucleoli Marker Antibody (NM95).




Immunohistochemistry-Paraffin: Nucleoli Marker Antibody (NM95) - Azide and BSA Free [NBP2-34695] - Formalin-fixed, paraffin-embedded human Tonsil stained with Nucleoli Marker Antibody (NM95).



Flow Cytometry: Nucleoli Marker Antibody (NM95) - Azide and BSA Free [NBP2-34695] - Flow Cytometric Analysis of PFA-fixed K562 cells using Nucleoli Marker Antibody (NM95) followed by Goat anti-Mouse IgG-CF488 (Blue); Isotype Control (Red)



Immunohistochemistry-Paraffin: Mouse Monoclonal Nucleoli Marker Antibody (NM95) - Azide and BSA Free [NBP2-34695] - Nucleoli Marker Antibody (NM95) [Biotin] (Catalog # NBP2-34695B) immunoreactivity in adipose surrounding mouse kidney. The antibody was diluted 1 to 100 and left on tissue sections for 16h at 4 degrees Celsius. Image from a verified customer review.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-34695-0.1mg

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-34695

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

