

# Product Datasheet

## Cytokeratin 3/2p Antibody (2Q1040(AE5)) [Biotin] NBP2-34844B

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-34844B](http://www.novusbio.com/NBP2-34844B)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-34844B](http://www.novusbio.com/reviews/destination/NBP2-34844B)



**NBP2-34844B**

Cytokeratin 3/2p Antibody (2Q1040(AE5)) [Biotin]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                     |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                      |
| Clone                       | 2Q1040(AE5)                                                                                                                                                                                                                                                                                                                                     |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                              |
| Isotype                     | IgG1                                                                                                                                                                                                                                                                                                                                            |
| Conjugate                   | Biotin                                                                                                                                                                                                                                                                                                                                          |
| Purity                      | Ion exchange chromatography                                                                                                                                                                                                                                                                                                                     |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                                             |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                 |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                     | 3850                                                                                                                                                                                                                                                                                                                                            |
| Gene Symbol                 | KRT3                                                                                                                                                                                                                                                                                                                                            |
| Species                     | Human, Mouse, Bovine, Rabbit                                                                                                                                                                                                                                                                                                                    |
| Reactivity Notes            | Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Specificity/Sensitivity     | Recognizes the 64kD polypeptide/keratin K3 (cytokeratin 3) of rabbit corneal epithelium and keratin K76 (cytokeratin K2p) of palate epithelium. As such represents an excellent marker for studying corneal epithelial-type differentiation. Also recognizes bovine snout and lip epithelium, as well as a rabbit corneal cell line.            |
| Immunogen                   | Rabbit corneal epithelial keratin                                                                                                                                                                                                                                                                                                               |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                 |
| Applications                | Western Blot, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                  |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen                                                                                                                                                                                                                                                  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                          |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP2-34844B**

---

|            |                                           |
|------------|-------------------------------------------|
| NBP1-96888 | Mouse IgG1 Isotype Control (MG1) [Biotin] |
| 236-EG-200 | EGF [Unconjugated]                        |
| NB100-355  | RPE65 Antibody (401.8B11.3D9) - BSA Free  |
| NBP2-29429 | Cytokeratin, pan Antibody (AE-1/AE-3)     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-34844B](http://www.novusbio.com/reviews/submit/NBP2-34844B)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

