

Product Datasheet

Recombinant Rat IL-1ra/IL-1F3/IL1RN Protein NBP2-35225-100ug

Unit Size: 100 ug

Store at -20 to -70C as supplied. After reconstitution, store at 2 to 8C for 1 month and at -20 to -70C for long term storage. Avoid repeated freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-35225

Updated 1/25/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-35225




NBP2-35225-100ug

Recombinant Rat IL-1ra/IL-1F3/IL1RN Protein

Product Information	
Unit Size	100 ug
Concentration	Lyoph
Storage	Store at -20 to -70C as supplied. After reconstitution, store at 2 to 8C for 1 month and at -20 to -70C for long term storage. Avoid repeated freeze-thaw cycles.
Preservative	No Preservative
Reconstitution Instructions	Recommended to centrifuge prior to opening. Reconstitute in sterile distilled water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0mg/mL.
Purity	>95%, by SDS-PAGE and HPLC
Buffer	Lyophilized from a 0.2 um filtered concentrated solution in PBS, pH 7.4.
Target Molecular Weight	17.4 kDa
Product Description	
Description	A single non-glycosylated polypeptide chain containing 152 amino acids corresponding to IL-1ra/IL-1F3/IL1RN Source: <i>E. coli</i> Uniprot ID: P25086 Amino Acid Sequence: HPAGKRPCCKM QAFRIWDTNQ KTFYLRNNQL IAGYLQGPNT KLEEKIDMVP IDFRNVFLGI HGGKLCCLSCV KSGDDTKLQL EEVNITDLNK NKEEDKRFTF IRSETGPTTS FESLACPGWF LCTTLEADHP VSLTNTPKPEP CTVTKFYFQE DQ
Gene ID	3557
Gene Symbol	IL1RN
Species	Rat
Details of Functionality	IL-1ra/IL-1F3/IL1RN protein is fully biologically active when compared to standard. The ED50 as determined by inhibiting IL-1 alpha-dependent proliferation of murine D10S cells is less than 150 ng/ml, corresponding to a specific activity of > 6.7 x 10 ³ IU/mg in the presence of 50 pg/ml IL-1ra/IL-1F3/IL1RN
Endotoxin Note	Less than 1 EU/ug of IL-1ra/IL-1F3/IL1RN as determined by LAL method.
Product Application Details	
Applications	SDS-Page, Bioactivity
Recommended Dilutions	SDS-Page, Bioactivity

Images

SDS-Page: Rat IL-1ra/IL-1F3/IL1RN Protein [NBP2-35225]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-35225-100ug

NBP2-35105-100ug	Recombinant Mouse IL-1ra/IL-1F3/IL1RN Protein
210-TA-005	TNF-alpha [Unconjugated]
DRA00B	IL-1ra/IL-1F3/IL1RN [HRP]
M6000B-1	IL-6 [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-35225

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

