

# Product Datasheet

## Histone H3 Antibody (1B1-B2) [Allophycocyanin] NBP2-36468APC

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-36468APC](http://www.novusbio.com/NBP2-36468APC)

Updated 10/29/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-36468APC](http://www.novusbio.com/reviews/destination/NBP2-36468APC)



**NBP2-36468APC**

Histone H3 Antibody (1B1-B2) [Allophycocyanin]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 1B1-B2                                                                                      |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG3 Kappa                                                                                  |
| Conjugate           | Allophycocyanin                                                                             |
| Purity              | Protein G purified                                                                          |
| Buffer              | PBS                                                                                         |

| Product Description |                                                                                                                             |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Host                | Mouse                                                                                                                       |
| Gene ID             | 126961                                                                                                                      |
| Gene Symbol         | H3C14                                                                                                                       |
| Species             | Human, Mouse, Rat                                                                                                           |
| Immunogen           | This Histone H3 antibody (1B1-B2) was raised against synthetic peptide corresponding to the C-terminus of human Histone H3. |

| Product Application Details |                                                                                                                                                                      |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Chromatin Immunoprecipitation (ChIP), Single Cell Western |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin, Chromatin Immunoprecipitation (ChIP), Single Cell Western |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-36468APC**

---

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| NB21-1021PEP       | Histone H3 [Monomethyl Lys4] Antibody Blocking Peptide  |
| H00008336-P01-10ug | Recombinant Human HIST1H2AM GST (N-Term) Protein        |
| NB21-1251PEP       | Histone H3 [Monomethyl Lys36] Antibody Blocking Peptide |
| NB500-171          | Histone H3 Antibody - BSA Free                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-36468APC](http://www.novusbio.com/reviews/submit/NBP2-36468APC)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

