

Product Datasheet

MTHFR Antibody (5D3) - BSA Free NBP2-37607

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-37607

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-37607



NBP2-37607

MTHFR Antibody (5D3) - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	5D3
Preservative	0.03% Sodium Azide
Isotype	IgG1
Purity	Unpurified
Buffer	Ascites
Target Molecular Weight	75 kDa

Product Description


Host	Mouse
Gene ID	4524
Gene Symbol	MTHFR
Species	Human, Rat
Immunogen	Purified recombinant fragment of human MTHFR expressed in E. Coli.

Product Application Details


Applications	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500 - 1:2000, ELISA 1:10000, Immunohistochemistry 1:200 - 1:1000, Immunohistochemistry-Paraffin

Images


Western Blot: MTHFR Antibody (5D3) [NBP2-37607] - Western blot analysis using MTHFR mouse mAb against Rat Heart cell lysate.




Immunohistochemistry-Paraffin: MTHFR Antibody (5D3) [NBP2-37607] - Analysis of intima cancer tissues (left) and prostate tissues (right) using MTHFR mouse mAb with DAB staining.




Western Blot: MTHFR Antibody (5D3) [NBP2-37607] - Western blot analysis using MTHFR mAb against HEK293 (1) and MTHFR(AA: 339-499)-hlgGfc transfected HEK293 (2) cell lysate.



Immunohistochemistry-Paraffin: MTHFR Antibody (5D3) [NBP2-37607] - Analysis of lung cancer using MTHFR mouse mAb with DAB staining.



ELISA: MTHFR Antibody (5D3) [NBP2-37607] - Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-37607

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP3-21401PEP	MTHFR Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-37607

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

