

Product Datasheet

CD2 Antibody (TS1/8) - BSA Free NBP2-37714

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-37714

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-37714



NBP2-37714

CD2 Antibody (TS1/8) - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	TS1/8
Preservative	15mM Sodium Azide
Isotype	IgG1 Kappa
Purity	Protein A purified
Buffer	Phosphate buffered saline (PBS), pH 7.4
Target Molecular Weight	50 kDa

Product Description

Description	Novus Biologicals Mouse CD2 Antibody (TS1/8) - BSA Free (NBP2-37714) is a monoclonal antibody validated for use in Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	914
Gene Symbol	CD2
Species	Human
Specificity/Sensitivity	The mouse monoclonal antibody TS1/8 recognizes CD2, a 50 kDa glycoprotein present on the human peripheral blood T lymphocytes and NK cells; also expressed by all thymocytes. HLDA V; WS Code S025
Immunogen	Cytotoxic T lymphocytes.


Product Application Details

Applications	Flow Cytometry, CyTOF-ready
Recommended Dilutions	Flow Cytometry 1-4 ug/ml, CyTOF-ready
Application Notes	This antibody is Cytof ready.




Images


Flow Cytometry: CD2 Antibody (TS1/8) [NBP2-37714] - Staining of CD2 in human peripheral blood with anti-CD2 (TS1/8) PE.




Flow Cytometry: CD2 Antibody (TS1/8) [NBP2-37714] - Staining of CD2 in human peripheral blood with anti-CD2 (TS1/8) PE-CyTM7.



Flow Cytometry: CD2 Antibody (TS1/8) [NBP2-37714] - Staining of CD2 in human peripheral blood with anti-CD2 (TS1/8) APC-CyTM7.



Flow Cytometry: CD2 Antibody (TS1/8) - BSA Free [NBP2-37714] - Surface staining pattern of human peripheral blood stained using anti-human CD2 (TS1/8) purified antibody (low endotoxin, concentration in sample 4 ug/ml) GAM APC.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-37714

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
NBP2-22745	Recombinant Human CD2 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-37714

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

