

Product Datasheet

TNF RI/TNFRSF1A Antibody (H398) [Allophycocyanin] NBP2-37724APC

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-37724APC

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-37724APC



NBP2-37724APC

TNF RI/TNFRSF1A Antibody (H398) [Allophycocyanin]

| Product Information | |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | H398 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a |
| Conjugate | Allophycocyanin |
| Purity | Protein A purified |
| Buffer | PBS |
| Product Description | |
| Host | Mouse |
| Gene ID | 7132 |
| Gene Symbol | TNFRSF1A |
| Species | Human |
| Specificity/Sensitivity | Antibody recognizes the extracellular domain of CD120a, a 55 kDa receptor for tumor necrosis factor. Antibody blocks biological activity of both natural and recombinant human TNF alpha and TNF beta. |
| Immunogen | Recombinant full length human CD120a. |
| Endotoxin Note | Endotoxin level is less than 0.01 EU/ug of the protein, as determined by the LAL test. |
| Product Application Details | |
| Applications | Flow Cytometry, Functional, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin, Immunoprecipitation, Block/Neutralize, CyTOF-ready |
| Recommended Dilutions | Flow Cytometry, Immunohistochemistry, Immunoprecipitation, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, Functional, CyTOF-ready, Block/Neutralize |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-37724APC

| | |
|----------------|-------------------------------------------------------|
| IC003A | Mouse IgG2a Isotype Control (20102) [Allophycocyanin] |
| NBP2-35274-5ug | Recombinant Human TNF RI/TNFRSF1A Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| DRT100 | TNF RI/TNFRSF1A [HRP] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-37724APC

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

