

Product Datasheet

V5 Epitope Tag Antibody (E10/V4RR) - Azide and BSA Free NBP2-37825

Unit Size: 100ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-37825

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-37825



NBP2-37825

V5 Epitope Tag Antibody (E10/V4RR) - Azide and BSA Free

Product Information	
Unit Size	100ug
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	E10/V4RR
Preservative	No Preservative
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2)
Product Description	
Host	Mouse
Species	Epitope Tag
Specificity/Sensitivity	This product recognizes V5 (GKPIPPELLGLDST) tagged fusion proteins. It has successfully been used in immunoprecipitation, ELISA, Western blotting, and immunofluorescence applications.
Immunogen	GKPIPPELLGLDST (V5) synthetic peptide conjugated to KLH
Product Application Details	
Applications	Western Blot, ELISA, Immunoprecipitation, Immunofluorescence
Recommended Dilutions	Western Blot 1:1000 - 1:3000, ELISA 1:100 - 1:2000, Immunoprecipitation 1:10 - 1:500, Immunofluorescence 1:500 - 1:2000




Images

Western Blot: V5 Epitope Tag Antibody (E10/V4RR) [NBP2-37825] - Analysis of 15ug of 293T cell lysates, from cells transfected with a V5-Wnt3a or V5-Vstm2 construct for the indicated number of hours, per well. Data courtesy of the Innovators Program.



Western Blot: V5 Epitope Tag Antibody (E10/V4RR) [NBP2-37825] - Analysis of 1:1000 (1ug/mL) Ab dilution probed against HEK293 cells transfected with V5-tagged protein vector; untransfected (1) and transfected (2).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-37825

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-37825

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

