

Product Datasheet

dynein, cytoplasmic 2, light intermediate chain 1 Recombinant Protein Antigen NBP2-38008PEP

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-38008PEP

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-38008PEP



NBP2-38008PEP

dynein, cytoplasmic 2, light intermediate chain 1 Recombinant Protein Antigen

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Preservative	No Preservative
Purity	>80% by SDS-PAGE and Coomassie blue staining
Buffer	PBS and 1M Urea, pH 7.4.
Target Molecular Weight	27 kDa

Product Description	
Description	<p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human DYNC2LI1.</p> <p>Source: <i>E. coli</i></p> <p>Amino Acid Sequence: PENDIGKLHAHSPMELWKKVYEKLFPPKSINTLKDIDKDPARDPQYAENEVDEMR IQKDLELEQYKRSSSKSWKQIEL</p> <p>Fusion Tag: N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G)</p> <p>This product is intended to be used as a blocking antigen for antibody competition assays. Any other use of this antigen is done at the risk of the user. The use of this product for commercial production is strictly prohibited. Please contact technical support if you have any questions.</p>
Gene ID	51626
Gene Symbol	DYNC2LI1
Species	Human

Product Application Details	
Applications	Antibody Competition
Recommended Dilutions	Antibody Competition 10 - 100 molar excess

Application Notes

This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP2-38008.

It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml.

For current lot information, including availability, please contact our technical support team click nb-technical@bio-techne.com

For further blocking peptide related information and a protocol, click [here](#).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-38008PEP

H00051626-P01-10ug	Recombinant Human dynein, cytoplasmic 2, light intermediate chain 1 GST (N-Term) Protein
AF2400	HNF-3 beta/FoxA2 Antibody [Unconjugated]
NBP2-38008	dynein, cytoplasmic 2, light intermediate chain 1 Antibody
H00004218-M02	RAB8A Antibody (3G1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-38008PEP

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

