

Product Datasheet

STAT3 Antibody - BSA Free

NBP2-38857

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-38857

Updated 2/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-38857



NBP2-38857

STAT3 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Host	Rabbit
Gene ID	6774
Gene Symbol	STAT3
Species	Human
Immunogen	This antibody was developed against a recombinant protein corresponding to amino acids: EKPMEIARIVARCLWEESRLLQTAATAAQGGQANHPTAAVVTEKQQMLEQHL QDVRKRVQDLEQKMKVVENLQDDDFDFNYKTLKSQGDMDLNGNNSVTRQK MQQLEQMLTALDQMRRSIVSELAGLLSAMEYVQK


Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunohistochemistry-Paraffin 1:50 - 1:200
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Western Blot: STAT3 Antibody [NBP2-38857] - Analysis in human cell lines A-549 and PC-3 using anti-STAT3 antibody. Corresponding STAT3 RNA-seq data are presented for the same cell lines. Loading control: anti-PPIB.




Immunohistochemistry-Paraffin: STAT3 Antibody [NBP2-38857] - Staining of human spleen shows moderate cytoplasmic positivity in cells in red pulp.




Immunohistochemistry-Paraffin: STAT3 Antibody [NBP2-38857] - Staining of human cervix, uterine shows moderate nuclear positivity in squamous epithelial cells.




Immunohistochemistry-Paraffin: STAT3 Antibody [NBP2-38857] - Staining of human fallopian tube shows moderate cytoplasmic and nuclear positivity in glandular cells.



Immunohistochemistry-Paraffin: STAT3 Antibody [NBP2-38857] - Staining of human skeletal muscle shows weak nuclear positivity in myocytes.



Immunohistochemistry-Paraffin: STAT3 Antibody [NBP2-38857] - Staining of human tonsil shows strong cytoplasmic positivity in squamous epithelial cells.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-38857

NBP2-38857PEP	STAT3 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-38857

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

