

# Product Datasheet

## Slain2 Antibody [PerCP] NBP2-41117PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-41117PCP](http://www.novusbio.com/NBP2-41117PCP)

Updated 10/24/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-41117PCP](http://www.novusbio.com/reviews/destination/NBP2-41117PCP)



**NBP2-41117PCP**

Slain2 Antibody [PerCP]

| Product Information         |                                                                                                                                                                                                            |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                     |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                   |
| Clonality                   | Polyclonal                                                                                                                                                                                                 |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                         |
| Isotype                     | IgG                                                                                                                                                                                                        |
| Conjugate                   | PerCP                                                                                                                                                                                                      |
| Purity                      | Peptide affinity purified                                                                                                                                                                                  |
| Buffer                      | PBS                                                                                                                                                                                                        |
| Product Description         |                                                                                                                                                                                                            |
| Host                        | Rabbit                                                                                                                                                                                                     |
| Gene ID                     | 57606                                                                                                                                                                                                      |
| Gene Symbol                 | SLAIN2                                                                                                                                                                                                     |
| Species                     | Human                                                                                                                                                                                                      |
| Specificity/Sensitivity     | SLAIN2 antibody is predicted not to cross-react with SLAIN1.                                                                                                                                               |
| Immunogen                   | Antibody was raised against a 14 amino acid synthetic peptide near the amino terminus of human SLAIN2. The immunogen is located within amino acids 70 - 120 of SLAIN2. Amino Acid Sequence: PRRTSSEELRDATS |
| Product Application Details |                                                                                                                                                                                                            |
| Applications                | ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                 |
| Recommended Dilutions       | ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                 |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                     |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

**Products Related to NBP2-41117PCP**

---

|                    |                                               |
|--------------------|-----------------------------------------------|
| NBP2-24891PCP      | Rabbit IgG Isotype Control [PerCP]            |
| H00057606-P01-10ug | Recombinant Human Slain2 GST (N-Term) Protein |
| 3218-ND-025        | Nodal [Unconjugated]                          |
| NBP2-57345PEP      | Slain2 Recombinant Protein Antigen            |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-41117PCP](http://www.novusbio.com/reviews/submit/NBP2-41117PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

