

Product Datasheet

IL-28B/IFN-lambda 3 Antibody - BSA Free NBP2-41236

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-41236

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-41236



NBP2-41236

IL-28B/IFN-lambda 3 Antibody - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	16 kDa

Product Description

Description	Novus Biologicals Rabbit IL-28B/IFN-lambda 3 Antibody - BSA Free (NBP2-41236) is a polyclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	282617
Gene Symbol	IFNL3
Species	Human, Porcine
Reactivity Notes	Porcine reactivity reported from a verified customer review.
Specificity/Sensitivity	IL-28B antibody is human specific. IL-28B antibody may cross-react with IL-28A and IL-29.
Immunogen	Antibody was raised against a 17 amino acid peptide near the center of human IL-28B. The immunogen is located within amino acids 120 - 170 of IL-28B . Amino Acid Sequence: IQPQPTAGPRTRGRLH


Product Application Details

Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot 1-2 ug/ml, ELISA, Immunocytochemistry/ Immunofluorescence 5 ug/ml, Immunohistochemistry-Paraffin Validated for IHC-P from a verified customer review.




Images


Western Blot: IL-28B/IFN-lambda 3 Antibody [NBP2-41236] - Analysis of IL-28B in HeLa cell lysate with IL-28B antibody at (A) 1 and (B) 2 ug/ml.



Immunohistochemistry-Paraffin: Rabbit Polyclonal IL-28B/IFN-lambda 3 Antibody [NBP2-41236] - Detection of IFN- λ 3 by immunofluorescence staining in a hepatic tissue of a conventional pig, showing a low amount of IFN- λ 3 in the perinuclear, cytoplasmic region of hepatocytes frequently forming the hepatic cords. x300. DAPI stained nuclei. Image from a verified customer review.



Immunocytochemistry/ Immunofluorescence: IL-28B/IFN-lambda 3 Antibody - BSA Free [NBP2-41236] - Immunofluorescence of IL-28B/IFN-lambda 3 in HeLa cells with IL-28B/IFN-lambda 3 antibody at 5 ug/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-41236

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-35077-5ug	Recombinant Mouse IL-28B/IFN-lambda 3 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-41236

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

