

Product Datasheet

Vinculin Antibody - BSA Free

NBP2-41237

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-41237

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-41237



NBP2-41237

Vinculin Antibody - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	124 kDa

Product Description


Host	Rabbit
Gene ID	7414
Gene Symbol	VCL
Species	Human, Mouse, Rat
Specificity/Sensitivity	Vinculin antibody is human, mouse and rat specific. At least three isoforms of Vinculin are known to exist. Human Vinculin has 3 isoforms, including isoform 2 (1134aa, 124kD), isoform 1 (1066aa, 117kD) and isoform 3 (222aa, 25kD). Mouse Vinculin has one isoform (1066aa, 117kD) and Rat Vinculin also has one isoform (1066aa, 117kD). NBP2-41237 can detect human, mouse and rat.
Immunogen	Antibody was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of human Vinculin. The immunogen is located within amino acids 1060-1110 of Vinculin. Amino Acid Sequence: LGRTNISDEESEQATE

Product Application Details


Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1-2 ug/mL, ELISA, Immunohistochemistry 5 ug/mL, Immunocytochemistry/ Immunofluorescence 20 ug/mL, Immunohistochemistry-Paraffin 1:10-1:500

Images


Western Blot: Vinculin Antibody [NBP2-41237] - Western Blot Validation in Mouse 3T3/balb Cell Lysate. Loading: 15 ug of lysates per lane. Antibodies: Vinculin NBP2-41237 (Lane 1: 0.5 ug/mL and Lane 2: 1 ug/mL), 1h incubation at RT in 5% NFDN/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.




Immunohistochemistry: Vinculin Antibody - BSA Free [NBP2-41237] - Vinculin Antibody [NBP2-41237] - Rat small intestine tissue with Vinculin antibody at 20 ug/mL.




Western Blot: Vinculin Antibody [NBP2-41237] - Independent Antibody Validation (IAV) via Protein Expression Profile in Human, Mouse and Rat Cell Lines. Loading: 15 ug of lysates per lane. Antibodies: Vinculin NBP2-41237 (2 ug/mL), Vinculin from Company B (1 ug/mL), Vinculin NBP2-41274 (1 ug/mL), and beta-actin (1.5 ug/mL), 1h incubation at RT in 5% NFDN/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.




Western Blot: Vinculin Antibody [NBP2-41237] - Western Blot Validation in Human PC-3 Cell Lysate. Loading: 15 ug of lysates per lane. Antibodies: Vinculin NBP2-41237 (Lane 1: 0.5 ug/mL and Lane 2: 1 ug/mL), 1h incubation at RT in 5% NFDN/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.



Western Blot: Vinculin Antibody [NBP2-41237] - Western Blot Validation in Human Cell Lines. Loading: 15 ug of lysates per lane. Antibodies: Vinculin NBP2-41237 (2 ug/mL), 1h incubation at RT in 5% NFDN/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.



Immunohistochemistry: Vinculin Antibody [NBP2-41237] - Rat small intestine tissue with Vinculin antibody at 5 ug/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-41237

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00007414-P01-10ug	Recombinant Human Vinculin GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-41237

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

