

# Product Datasheet

## GDPD5 Antibody [PerCP] NBP2-41268PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-41268PCP](http://www.novusbio.com/NBP2-41268PCP)

Updated 10/24/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-41268PCP](http://www.novusbio.com/reviews/destination/NBP2-41268PCP)



**NBP2-41268PCP**

GDPD5 Antibody [PerCP]

| <b>Product Information</b>         |                                                                                                                                                                                                       |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                                |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                           |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                              |
| <b>Clonality</b>                   | Polyclonal                                                                                                                                                                                            |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                    |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                   |
| <b>Conjugate</b>                   | PerCP                                                                                                                                                                                                 |
| <b>Purity</b>                      | Peptide affinity purified                                                                                                                                                                             |
| <b>Buffer</b>                      | PBS                                                                                                                                                                                                   |
| <b>Product Description</b>         |                                                                                                                                                                                                       |
| <b>Host</b>                        | Rabbit                                                                                                                                                                                                |
| <b>Gene ID</b>                     | 81544                                                                                                                                                                                                 |
| <b>Gene Symbol</b>                 | GDPD5                                                                                                                                                                                                 |
| <b>Species</b>                     | Human, Mouse, Rat                                                                                                                                                                                     |
| <b>Specificity/Sensitivity</b>     | At least three isoforms of GDPD5 are known to exist.                                                                                                                                                  |
| <b>Immunogen</b>                   | Antibody was raised against an 18 amino acid peptide near the amino terminus of human GDPD5. The immunogen is located within the first 50 amino acids of GDPD5. Amino Acid Squence: RWKRYQRSHDDTTPWER |
| <b>Product Application Details</b> |                                                                                                                                                                                                       |
| <b>Applications</b>                | Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry                                                                                      |
| <b>Recommended Dilutions</b>       | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                                                                                      |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-41268PCP**

---

|               |                                    |
|---------------|------------------------------------|
| NBP2-24891PCP | Rabbit IgG Isotype Control [PerCP] |
| NBP2-56359PEP | GDPD5 Recombinant Protein Antigen  |
| DY8198-05     | alpha-Defensin 1 [Biotin]          |
| NBP2-49438PEP | GDPD5 Recombinant Protein Antigen  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-41268PCP](http://www.novusbio.com/reviews/submit/NBP2-41268PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

