

Product Datasheet

HIV-1 Gag p24 Antibody (8G9) - BSA Free NBP2-41336

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-41336

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-41336



NBP2-41336


HIV-1 Gag p24 Antibody (8G9) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	8G9
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS
Target Molecular Weight	24, 41, 55 kDa
Product Description	
Host	Mouse
Gene ID	155030
Gene Symbol	gag
Species	Virus
Specificity/Sensitivity	By Western blot, anti-HIV-1 Gag p24 antibody detects a ~24 kDa, a ~41 kDa, and a ~55 kDa protein, corresponding to HIV-1 Gag p24 and to its precursors p41 and p55, respectively, in HIV-1 samples.
Immunogen	Antibody was raised against a recombinant full-length HIV-1 Gag p24 protein. Amino Acid Sequence: pivqniqggqmvhqaisprtlnawvkvveekafspevipmfsalsegatpqdlntmlntvgghqaamqmlketi neeeaaewdrvhpvhagpiapgqmreprgsdiagttstlqeqigwmtnnppipvgeiykrwiilglnkivrmys ptsildirqgpkpfrdyvdrfyktraeqasqevknwmtetllvqnanpdcktilkalgpaatleemmtacqgv gpghkarvla
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 0.2-0.5 ug/ml, ELISA 0.2-0.5 ug/ml




Images

Western Blot: HIV-1 Gag p24 Antibody (8G9) - BSA Free [NBP2-41336]
- Western blot analysis of 20 ng of (A) viral p24 and (B) recombinant p24 with anti-HIV-1 Gag p24 antibody NBP2-41336 at (A) 0.5 ug/mL and (B) 0.2 ug/mL, respectively.



ELISA: HIV-1 Gag p24 Antibody (8G9) - BSA Free [NBP2-41336] -
Titration ELISA of anti-HIV-1 Gag p24 antibody NBP2-41336 with 100 ng of recombinant HIV-1 Gag p24 protein.



Publications

Puertas MC, Bayon-Gil A, Garcia-Guerrero MC Et al. VIP-SPOT: an Innovative Assay To Quantify the Productive HIV-1 Reservoir in the Monitoring of Cure Strategies mBio 2021-06-22 [PMID: 34154408]

Details:

This citation used the Biotin format of this antibody.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-41336

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP1-46036-0.1mg	Recombinant HIV-1 Gag p24 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-41336

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

