

# Product Datasheet

## Hsp70 interacting protein HIP Antibody (5C6) [DyLight 650] NBP2-42671C

Unit Size: 100 ul

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-42671C](http://www.novusbio.com/NBP2-42671C)

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-42671C](http://www.novusbio.com/reviews/destination/NBP2-42671C)



**NBP2-42671C**

Hsp70 interacting protein HIP Antibody (5C6) [DyLight 650]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 100 ul                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 5C6                                                                                         |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG2a Kappa                                                                                 |
| Conjugate           | DyLight 650                                                                                 |
| Purity              | Protein A purified                                                                          |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description |                                                                                                                                                                                       |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |
| Host                | Mouse                                                                                                                                                                                 |
| Gene ID             | 6767                                                                                                                                                                                  |
| Gene Symbol         | ST13                                                                                                                                                                                  |
| Species             | Human                                                                                                                                                                                 |
| Immunogen           | Recombinant human Hsp70 interacting protein HIP (1-369aa) purified from E. coli                                                                                                       |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                     |

| Product Application Details |                                                                             |
|-----------------------------|-----------------------------------------------------------------------------|
| Applications                | Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence |
| Recommended Dilutions       | Western Blot, Flow Cytometry, ELISA, Immunocytochemistry/Immunofluorescence |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-42671C**

---

|               |                                                           |
|---------------|-----------------------------------------------------------|
| NBP1-96981C   | Mouse IgG2a Kappa Isotype Control (M2AK) [DyLight 650]    |
| NBC1-19066    | Recombinant Human Hsp70 interacting protein HIP Protein   |
| NBP2-47427PEP | Hsp70 interacting protein HIP Recombinant Protein Antigen |
| NB300-164     | Bestrophin 1 Antibody (E6-6) - BSA Free                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-42671C](http://www.novusbio.com/reviews/submit/NBP2-42671C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

